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There is no defined protocol that we recommend for tissue extraction as this may differ for different tissue types and/or analyte one may wish to measure. Generally, most protocols that are used in ELISAs can be used, but here are some guidelines in selecting a method.

- 1) Homogenize cells or tissues mechanically (eg. Ultrasonication) in a PBS based buffer containing protease inhibitors (like aprotinin or an inhibitor cocktail) and low (< 0.2%) non-ionic detergent concentration.
- 2) Extraction medium should not contain any organic solvents like DMSO etc.
- 3) Centrifuge the extract and store the supernatant at < -20 C.
- 4) Use the extraction medium as matrix in blank, standard curve and QCs.

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1. [Insulinitis and characterisation of infiltrating T cells in surgical pancreatic tail resections from patients at onset of type 1 diabetes](#). Lars Krogvold et al., *Diabetologia*. 2016; 59:492–501.
2. [Nilotinib and bosutinib modulate pre-plaque alterations of blood immune markers and neuro-inflammation in Alzheimer's disease models](#). I. Lonskaya et al., *Neuroscience*. 2015 Sep 24;304:316-27.
3. [Cranial irradiation significantly reduces beta amyloid plaques in the brain and improves cognition in a murine model of Alzheimer's Disease \(AD\)](#). Brian Marples et al., *Radiother Oncol*. 2016; 118(3):579-80.
4. [Helicobacter infection alters MyD88 and Trif signalling in response to intestinal ischaemia-reperfusion](#). Hoffman SM et al., *Exp Physiol*. 2011; 96(2):104-13.
5. [Nod-like receptor pyrin containing 3 \(NLRP3\) in the post-mortem frontal cortex from patients with bipolar disorder: A potential mediator between mitochondria and immune-activation](#). HK Kim et al., *Journal of Psychiatric Research*, Volume 72, January 2016, Pages 43-50
6. [Aldose reductase regulates doxorubicin-induced immune and inflammatory responses by activating mitochondrial biogenesis](#). H Sonowal et al., *European Journal of Pharmacology*, Volume 895, 15 March 2021, 173884
7. [The JAK2 Inhibitor AZD1480 Potently Blocks Stat3 Signaling and Oncogenesis in Solid Tumors](#). Michael Hedvat et al., *Cancer Cell*, Volume 16, Issue 6, 8 December 2009, Pages 487-497

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本紙記載の製品は試験・研究用です。ヒト、動物への治療、もしくは診断目的として使用しないようご注意ください。なお、品目、製品情報、価格等は予告なく変更される場合がございます。予めご了承ください。記載内容は2023年2月時点の情報です。©2023 Merck KGaA, Darmstadt, Germany. All rights reserved. Merck, the vibrant M, and Millipore are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.