

# Enveloped Virus-Like Particle (VLP) Process Development

R. Silva<sup>1,2</sup>, R. Castro<sup>1,2</sup>, M.Q. Sousa<sup>1,2</sup>, H.R. Soares<sup>1,2</sup>, A.S. Coroadinha<sup>1,2</sup>, A. Stein<sup>3</sup>, S. Ribaud<sup>4</sup>, A. Xenopoulos<sup>5</sup>, P.M. Alves<sup>1,2</sup>, C. Peixoto<sup>1,2</sup>, M.J.T. Carrondo<sup>1,2,6</sup>

<sup>1</sup>IBET-Instituto de Biologia Experimental e Tecnológica, Oeiras, Portugal. <sup>2</sup>Instituto de Tecnologia Química e Biológica – UNL, Oeiras, Portugal. <sup>3</sup>Merck KGaA, Frankfurter Str. 250 D 64293 Darmstadt, Germany.

<sup>4</sup>Merck Millipore, 39 Route Industrielle de la Hardt, 67120 Molsheim, France. <sup>5</sup>EMD Millipore, 80 Ashby Road, Bedford, MA 01730, USA. <sup>6</sup>FCT/UNL, P-2829-516 Caparica, Portugal.

## Introduction

Hepatitis C virus (HCV) infection is a major public health problem, causing more than 350.000 deaths every year, according to the World Health Organization. Currently, there is no Hepatitis C vaccine and the standard treatment for acute HCV infections has several limitations, including its low efficacy<sup>1,2</sup>. There is thus a clear need for the development of a vaccine with both preventive and therapeutic roles.

Virus-like particles (VLPs) have received increased attention as vaccine candidates, since some of those already marketed have been very successful.

VLP titers have been helped by advances in systems biotechnology, as well as the flexibility and versatility of retroviral-like particles to incorporate and display proteins. In this work, we highlight these advances for a chimeric retro VLP displaying E1/E2 HCV proteins, produced in an insect cell expression system. This VLP is a candidate vaccine for Hepatitis C. Furthermore, rational tools are being developed to understand the purification processes of enveloped VLPs<sup>3</sup>.


In this study, the production and purification of VLP-HCV was investigated and improved. In particular disposable bioreactor technology was evaluated and compared with glass stirred tank bioreactor culture in terms of viable cell concentration, % viability, growth kinetics and stability.

Cell culture supernatants from both processes were then harvested and purified to understand how production affects the downstream process and product quality. Downstream train was improved by selecting an appropriate anion exchange (AEX) prototype resin yielding 70% recovery and a satisfactory baculovirus (BV) Log reduction (LRV) in bind and elute mode (B/E). Moreover, negative (flow-through, F-T) mode was also investigated, yielding a comparable BV clearance. Appropriate depth filtration and ultrafiltration (UF) devices were selected, leading to a GMP-compliant process that allows easy transfer to pilot scale level.

## Aim

Move towards the development of effective vaccine candidates for infectious diseases based on enveloped VLPs, using disposable and scalable technologies.

## Strategy



## Materials and Methods

**Cell lines and culture conditions.** *Spodoptera frugiperda* Sf9 cells (Life Technologies) were routinely cultivated in SF-900II SFM medium (Life Technologies) in shake-flasks at 27 °C and a stirring speed of 100 rpm.

**Determinations of hydrodynamic parameters.** Determination of the hydrodynamic parameters for the two bioreactors was performed as described by Cruz *et al.* 1998 and is shown in Table 1<sup>4</sup>.


**Bioreactor cultivations.** For VLP-HCV production in stirred tank bioreactor, Sf9 cells were cultured in 2 L working volume reusable stirred tank bioreactor (Biostat® B-DCU, Sartorius Stedim Biotech) or disposable Mobius® CellReady 3L Bioreactor (Merck Millipore). Cultivations were carried out at 27 °C, pO<sub>2</sub> of 30% and an aeration rate of 0.01 vvm. 24h after inoculation (1x10<sup>6</sup> cells/ml) cells were co-infected at an MOI of 2 for BV-Gag and BV-HCV, and culture supernatant was collected every day until bioreactor bulk harvesting at 96 hpi.

**Western-blot analysis.** VLP-HCV samples were separated on 4–12% acrylamide NuPAGE gels (Invitrogen). After transference, the PVDF membranes were incubated with anti-MLV-Gag (produced by hybridoma R187-ATCC CRL-1219), anti-HCV-E1 (Acris) or anti-HCV-E2 (Austral biologicals). The anti-mouse-HRP (GE Healthcare) or anti-RAT Alkaline Phosphatase (GE Healthcare) secondary antibodies were used and detection was performed using the ECL kit (GE Healthcare) or NBT/BCIP reagent (Thermo scientific), respectively.

**Clarification/Concentration/Purification of bioreactor bulk.** Bioreactor was clarified using different Optiscale® filter capsules with Polygard® CN filters (Merck Millipore). Concentration of the clarified bulks was performed with 50 cm<sup>2</sup> Pellicon® XL ultrafiltration cassettes (Merck Millipore). Six different ion-exchange prototypes (Merck Millipore) were evaluated.


## Upstream Processing of VLP-HCV

### I. Optimization of Sf9 growth in Mobius® CellReady 3L Bioreactor




### II. VLP-HCV production in bioreactor: reusable versus disposable vessels

#### Sf9 growth profile in the two bioreactors



#### VLP-HCV Productivity in the two bioreactors



## References

- [1] Manns, M.P., et al. (2006). Gut, 55(9): p. 1350-9.
- [2] Swadlow L, Klenerman P, Barnes E. Ever closer to a prophylactic vaccine for HCV. Expert Opin Biol Ther. Aug 2013; 13(8): 1109-1124.
- [3] Vicente, Tiago, et al. Analysis of adsorption of a baculovirus bioreaction bulk on an ion-exchange surface by surface plasmon resonance. J Biotech 148.4 (2010): 171-181
- [4] Cruz PE, Peixoto CC, Moreira JL, Carrondo MJ (1998). Biotechnol Bioeng, 60:408-18.

## Acknowledgements

The authors acknowledge the financial support received from Fundação para a Ciência e a Tecnologia, Portugal (PTDC/EBB-BIO/102649/2008 and EXPL/BBB-BIO/1129/2013). R. Castro and H.R. Soares acknowledge FCT for the award of the research grants SFRH/BPD/72523/2010 and SFRH/BD/81598/2011, respectively.




## Downstream Processing

### III. Clarification and Concentration Steps


Different clarification strategies with Optiscale® filter capsules with Polygard® CN filters were evaluated:

- Processes A to C use a combination of filters with different pore sizes;
- Process D couples bulk centrifugation to 0.6 µm device.

✓ Removal of 70% of DNA with strategy D.



### IV. Purification Step



- ✓ DBC<sub>10%</sub> of tested prototypes 1.1–1.3 µg HCV-VLP;
- ✓ Prototype A and B yielded a LRV reduction of 1.2 and a recovery of 65% and 70% respectively, for a B/E operation;

## Final Remarks

- Similar Sf9 growth profile was observed for the reusable STR and the disposable CR bioreactors;
- Comparable VLP-HCV production kinetics, as well as final productivity was obtained in the two types of bioreactors used;
- Development of a clarification and concentration train yielding 80% of recovery;
- Identification of the best IEX prototypes for B/E and F-T applications;

Merck Millipore is a division of

Merck Millipore, the M logo, Mobius, Optiscale, Polygard, and Pellicon are registered trademarks of Merck KGaA, Darmstadt, Germany. Lit Number: PS1923EN00 PS-14-10407 8/2014 © 2014 EMD Millipore Corporation, Billerica, MA, USA. All rights reserved.

Presented at

**ECI Vaccine Technology V Conference**  
Playa del Carmen, Mexico June 8–13, 2014