

2026-27

Analyte Quarterly

Reliable Biomarker Solutions

Be confident in your discoveries. From screening and qualified multiplexing to ultrasensitive biomarker detection and single protein analysis, all compatible with automation, you can rely on our platforms and services to accelerate your biomarker research.




MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

Millipore®

Preparation, Separation,
Filtration & Testing Products

Beyond Biomarkers

We are committed to delivering the best possible immunoassays for your research needs. Whether you are using our industry-leading MILLIPLEX® multiplex panels to broadly survey multiple analytes, performing femtogram/mL biomarker analysis using the FemtoQuest® ultrasensitive immunoassay platform, or examining single proteins with our tried-and-true ELISAs and RIAs, you can confidently expect reliable measurements, a simplified user experience, and knowledgeable scientific partners to help drive your biomarker project from hypothesis to publication.



Designed Reliability

Our Conferma® ELISAs are engineered for endogenous detection and lot-to-lot consistency of routine-use biomarkers to keep your data reproducible, run after run.



Genuine Flexability

Our industry-leading MILLIPLEX® multiplex assays for the Luminex® platform give you the FlexAbility to select your combination of biomarkers to simultaneously analyze with ease.




Next-Level Detectability

Our ultrasensitive SMC® kits for the FemtoQuest® platform take your research to the next level by detecting biomarkers down to fg/mL levels.



Complete Compareability

Our Belysa® Immunoassay Curve Fitting Software helps you quickly compare your standard and sample curves with ease, confirming that your methods and plates ran consistently.



Expert Customability

Our Custom Assay Development & Innovation (CADI) team of experienced scientists partner with you to develop custom immunoassays for your biomarker(s) of interest, accelerating your projects from discovery to clinical trials.



Keep your focus on the next stage.



Immunoassays at a Glance

		MILLIPLEX® Multiplex Assays		SMC® Assays	ELISAs
		Qualified Assays	Screening Assays		
Species	Human	•	•	•	•
	Mouse	•		•	•
	Rat	•		•	•
	Non-Human Primate	•		•	•
	Canine	•			•
	Feline	•			•
	Bovine	•			•
	Chicken	•			
	Equine	•			•
	Ovine	•			•
	Porcine	•			•
	Multi-Species	•		•	•
Analytes	Human	>620	115	>40	>1,000
	Mouse	>280		3	>300
	Rat	>170		3	>190
	Non-Human Primate	>65		1	>70
	Canine	>40		1	>30
	Feline	>15			>20
	Bovine	>15			>20
	Chicken	>10			
	Equine	>25			>20
	Ovine	>15			>5
	Porcine	>30			>30
	Multi-Species	>5		1	>5
Research Areas	Immunology/Inflammation	•	•	•	•
	Immune Response	•	•	•	•
	Virology	•		•	•
	Metabolism/Endocrinology	•	•	•	•
	Cardiovascular Disease	•	•	•	•
	Bone	•	•		•
	Cancer	•	•	•	•
	Neuroscience	•	•	•	•
	Toxicity	•	•		•
	Cell Signaling	•			•
Kit Format *only applies to MILLIPLEX® kits **Conferma® ELISAs only	Configurable (Select Your Own Analytes)*	•	•	-	-
	Premixed (Fixed Panel)*	•	•	-	-
	Bulk Premixed (Space Saver)*	•	•	-	-
	Custom Premix*	•	•	-	-
	Lot Matched at Kit Level			•	•**
Qualified Sample Types (see specific kits for more details)	Serum/Plasma	•		•	•
	Cell Culture Supernatants	•	Screening assays are sample agnostic	•	•
	Cell/Tissue Lysates	•			•
	Cerebrospinal Fluid	•		•	•
	Urine	•			•
Custom Kits and Services	—	•		•	•

New Immunoassay Products

New Instruments and Kits

Description	Species	Cat. No.
Instruments		
FemtoQuest® System	-	95-0200-00
SMC® High Sensitivity Kits		
2-Plex IL-17A, IL-22	Human	03-0225-2P
MILLIPLEX® Multiplex Panels		
Complement Panel	Rat/Mouse	RMCMP-260K
Cytokine/Chemokine/Growth Factor Expanded Panel 32-Plex	Porcine	PIGCYT-250K
Dual Reporter IgG/IgM Reagent Pack	-	DRAABPK-240K
Vector-Borne Disease Panel IgG	Human	HVBDG-270K

Coming Soon

Description	Species
SMC® High Sensitivity Kits	
pUB-S65	Human, Mouse, Rat
2-Plex Mutant HTT, Total HTT	Human
MILLIPLEX® Multiplex Panels	
Kidney Injury Panel	Non-Human Primate

Open Ultrasensitivity

Experience immunoassay flexibility with the ultrasensitive FemtoQuest® system

Measure down to fg/mL levels of biomarkers with verified SMC® technology and when your research needs require more, freely choose to chart your own course. Either customize your assay with our comprehensive all-in-one assay development kits or take it even further with our expert scientists for dual-plex custom assays, verification, and manufacturing.




Explore more at
SigmaAldrich.com/femtoquest



SMC® Technology

SMC® technology provides maximum immunoassay performance while following a workflow similar to traditional ELISA technology, as shown below for dual-plex analysis. By combining a unique assay elution step and robust digital counting on the FemtoQuest® system, SMC® technology improves signal-to-noise ratios over traditional immunoassay technologies, enhancing quantification at both low and high levels of expression with one system. The FemtoQuest® instrument can be used for measuring both single analyte and dual-plex assays.



Steps 2-3

Antibodies translate each biomarker into a signal

Step 4

Fluorescent dye-labeled detection antibodies dissociate from the immunocomplex

Step 5

On the FemtoQuest® instrument, signal is collected from free-floating fluorochromes using a rotating laser, camera, and counting individual molecules by an avalanche photodiode

Digital counting improves sensitivity and dynamic range

The FemtoQuest® system captures all digital events counted above the background threshold. Using enhanced algorithms, digital events are computed across a spectrum of time series for a single standard. This digital counting of single molecules improves the assay sensitivity and extends the dynamic range beyond what can be achieved with traditional methods.

Read More About Ultrasensitive Technology

Explore all the ways researchers are using SMC® ultrasensitive assays at SigmaAldrich.com/smc-learn



SMC® Technology for High Sensitivity Analysis

Advance Your Biomarker Quest



Achieve your desired immunoassay sensitivity with ease using the ultrasensitive FemtoQuest® system.

- Measure down to femtogram/mL levels with verified SMC® technology
- Easily perform dual-plex and single-analyte analysis with validated off-the-shelf kits
- Open platform to develop your own assays with all-in-one development kits
- Go fully custom with our CADI scientists
- Automation-ready and compatible with a variety of sample types
- Integrated data analysis with Belysa® software
- Low-maintenance, lightweight, fluidics-free benchtop instrument

Cat. No. 95-0200-00

Dynamic Range	>4 logs
Limit of Detection (LOD)	<1 fM
Instrument Footprint	Height: 358.53 mm (14.11 in) Width: 444.5 mm (17.5 in) Depth: 609.91 mm (24.0 in)
Instrument Weight	29.7 kg (65.5 lb)
Read Format	96-well assay format read on RFID-tagged SMC® Read Plates
Regulatory Compliance	LIMS compatible
Automation Compatible	Compatible with automation, including the AAW™ Workstation

Three Levels of Onsite Service Support for FemtoQuest® System

- **Total Service:** Unlimited visits including labor, travel, expenses and all parts for any service for one year.
- **Advanced Service:** Labor, travel and expenses for one preventative maintenance visit and one floating visit for one year.
- **Essential Service:** Labor for one preventative maintenance visit for one year.

Learn All About the Next-Gen Instrument

For more details about the FemtoQuest® system, visit SigmaAldrich.com/femtoquest-system



SMC® 2-Plex High Sensitivity Kits for Use on the FemtoQuest® Instrument

Product Description	Cat. No.	LLOQ (fg/mL) ¹	Species ²	Sample Type ³
Immunology/Inflammation				
2-Plex IL-1 β , IL-1 α	03-0218-2P	40 (IL-1 β), 40 (IL-1 α)	H	P,S
2-Plex IL-17A, IL-17F	03-0214-2P	20 (IL-17A), 90 (IL-17F)	H	P,S
2-Plex IL-17A, IL-22 NEW	03-0225-2P	35 (IL-17A), 87 (IL-22)	H	P,S
2-Plex IL-23, IL-22	03-0215-2P	40 (IL-23), 90 (IL-22)	H	P,S
2-Plex TNF α , IL-6	03-0216-2P	20 (TNF α), 60 (IL-6)	H	P,S
Cardiovascular				
2-Plex cTnI, IL-6	03-0217-2P	350 (cTnI), 50 (IL-6)	H	P,S

SMC® Singleplex High Sensitivity Kits

Product Description	Cat. No.	LLOQ (fg/mL) ¹	Species ²	Sample Type ³
Immunology/Inflammation				
IFN- α 2	03-0186-00	62	H	P,S
IFN- β 1	03-0191-00	150	H	P,S
IFN- γ	03-0181-00	33	H	P,S
IFN λ 1	03-0213-00	260	H	P,S
IL-1 α 2	03-0167-00	120	H	P,S
IL-1 β	03-0160-00	200	H	P,S
IL-2	03-0195-00	90	H	P,S
IL-4	03-0161-00	40	H	P,S
IL-5	03-0168-00	40	H	P,S
IL-6	03-0155-00	80	H	P,S
IL-10	03-0176-00	100	H	P,S
IL-12p70	03-0182-00	100	H	P,S
IL-13	03-0156-00	40	H	P,S
IL-17A/F Heterodimer	03-0177-00	2,340	H	P,S
IL-17A	03-0159-00	70	H	P,S
IL-17F	03-0164-00	200	H	P,S
IL-18	03-0190-00	90	H	P,S
IL-22	03-0162-00	100	H	P,S
IL-22BP	03-0219-00	3,910	H	P,S
IL-23	03-0157-00	100	H	P,S
MCP-1	03-0187-00	150	H	P,S
TNF α	03-0163-00	200	H	P,S
VEGF-A	03-0180-00	200	H	P,S

Get All the Best Practices

Our best practices article provides tips and tricks for using the FemtoQuest® system with SMC® assays, visit SigmaAldrich.com/femtoquest-best




SMC® Singleplex High Sensitivity Kits (continued)

Product Description	Cat. No.	LLOQ (fg/mL) ¹	Species ²	Sample Type ³
Neuroscience				
Phospho-α-Synuclein(Ser129)	03-0188-00	100	H	P,S,C
Total-α-Synuclein	03-0196-00	1,950	H	P,S,C
Amyloid beta 1-40	03-0145-00	5,860	H, M, R	P,C
Amyloid beta 1-42	03-0146-00	980	H, M, R	P,C
BDNF	03-0171-00	40	H	P,S,C
GFAP	03-0203-00	230	H	P,S,C
NF-L	03-0202-00	870	H	P,S,C
NPTX2	03-0199-00	580	H	P,S,C
SNAP-25	03-0206-00	260	H	P,S,C
TDP-43	03-0205-00	5,210	H	P,S,C
Phospho-Tau181	03-0184-00	940	H	P,S,C
Phospho-Tau217	03-0210-00	1,950	H	P,S,C
Phospho-Tau231	03-0211-00	1,950	H	P,S,C
Total Tau	03-0185-00	4,690	H	P,S,C
UCHL1	03-0183-00	2,600	H	P,S,C
Cancer				
CTLA-4	03-0220-00	300	H	P,S
PD-1	03-0207-00	3,470	H	P,S
PD-L1	03-0208-00	780	H	P,S
Cardiovascular				
cTnI	03-0154-00	690	H, Cy, R, C, GP	P,S
Metabolism/Endocrinology				
Glucagon	03-0153-00	780	H, M, R	P
Immunogenicity				
IgE ⁴	03-0201-00	31,000	H	P,S
IgM ⁴	03-0204-00	206,000	H	P,S
Virology				
SARS-CoV-2 RBD IgG ⁵	03-0193-00	NA ⁵	H	P,S
SARS-CoV-2 S1 IgA ⁵	03-0197-00	NA ⁵	H	P,S

¹ LLOQ: Lowest point on standard curve with CV <20% and accuracy within 20% of expected values. Lot-specific LLOQ will be reported on the Certificate of Analysis. ² Optimized for first species type listed. Other listed species have been tested, but not optimized for peak performance. Key: H=Human; M=Mouse; R= Rat; GP=Guinea Pig; Cy=Cynomolgus Monkey; C=Canine. ³ Optimized for use in sample type(s) listed. Key: P=Plasma; S=Serum; L=Lysate; U=Urine; C=Cerebrospinal Fluid (CSF). ⁴ Plate-based assay. ⁵ Qualitative assay.

What's In A Typical SMC® Kit?



Custom Assay Development Solutions

Advance your biomarker research from target to reliable insights.

Every lab works differently. With two flexible but compatible routes for creating an assay to measure your protein of interest, you can choose the approach that best suits your needs.

SMC® Assay Development Kits and Reagents for the FemtoQuest® Platform

Design custom SMC® immunoassays for your protein target and achieve femtogram-level sensitivity. When your analyte isn't available in an off-the-shelf kit, the FemtoQuest® platform plus companion kits and reagents enable you to build it quickly, flexibly, and with confidence.

Product Description	Cat. No.
SMC® Bead-Based Assay Development Kit	03-0222-00
SMC® Assay Buffer Optimization Kit	03-0221-00

Custom Assay Development & Innovation (CADI) Services

Advance your biomarker program without slowing down your science. The CADI team delivers custom SMC® and MILLIPLEX® immunoassays using a structured, scalable approach to ensure high-quality and reproducible results. From discovery to clinical research, we handle development so you can focus on what matters most.



From design through verification, our milestone-based approach with built-in data reviews at every stage keeps your project aligned, on schedule, and delivered to specification.



In addition to our assay development services, we can support your other custom assay needs:

Switching Formats

- MILLIPLEX® Assays: 96-well to 384-well
- SMC® Assays: bead-based to plate-based
- Converting MILLIPLEX® assays to SMC® assay platform

Sample Testing

- Prototype sample testing
- Custom kit sample testing (SMC® and MILLIPLEX® assays)
- Off-the-shelf immunoassay kit testing (SMC® and MILLIPLEX® assays)

Biotherapeutic Assay Development

- Pharmacokinetics
- Pharmacodynamics
- Immunogenicity

SMCxPlore™ Proof of Concept Sample Testing:

Use our SMCxPlore™ proof of concept sample testing service to “explore” SMC® technology with your analyte of interest or critical samples. Test up to 30 of your critical samples using one of our highly verified SMC® immunoassays.

From early screening through full-scale manufacturing, our comprehensive suite of tools and services supports every step of your biomarker journey.

Get started today!

View a complete listing of services at SigmaAldrich.com/customassay



Multiplexing Options That Fit Your Purpose

Introducing Screening and Qualified MILLIPLEX® Assays

Compelling PlexAbility with MILLIPLEX® Screening Assays


- Save sample, time, and money by profiling up to 115 analytes simultaneously without sacrificing specificity
- Excludes Serum Matrix to create sample agnostic panel for greater flexibility
- Ideal for detecting differences between experimental groups
- Human assay available
- Uses same critical raw materials as Qualified assays

Screen to identify analytes of interest and move to Qualified for further analysis



Genuine PlexAbility with MILLIPLEX® Qualified Assays

- Lot-to-lot reproducibility with rigorous QC testing
- Includes key reagents like Quality Controls and Serum Matrix for sample normalization
- Ideal for long-term or cross-site studies
- Broadest selection of analytes across research areas and species
- Configurable, premixed, and bulk kit options to fit your exact research needs



MILLIPLEX® Analyte Index

Analyte	Page No.	Analyte	Page No.	Analyte	Page No.
4-6					
4-1BBL/TNFSF9	21,25	Apolipoprotein E (Apo E)	21,22	Cardiac Troponin T (cTnT)	21,25,26
4E-BP1 (total)	30	ApoE4	22	Cathepsin D	21,22
4E-BP1 (Thr37/46)	30	APRIL	18,21	Caveolin-1 (CAV-1)	26
5'-NT/CD73	21,23,25,26	AQP4 (autoantibody)	22	CCHFV gp38	20
6Ckine/CCL21/Exodus-2	17,18	Arginase-1 (ARG1)	21,23,26	CCL25/TECK	23
A		Argonaute-2 (Ago2)	21	CCL28	18,19
α 1-Acid Glycoprotein (AGP)	21,25,26	Argonaute-2 (Ago2) (autoantibody)	22	CCL3L1	27
α-1-Antitrypsin (A1AT)	22	ATF2 (Thr71)	30	CCNB1/Cyclin B1 (autoantibody)	21
α-1 Microglobulin	23	ATR (total)	30	CD3ε (panTyr)	31
α-2 Macroglobulin (A2M)	21,22,25,26	AXL (soluble)	22	CD14 (soluble)	21
α-Fetoprotein (AFP)	20,21	B		CD25/IL-2Rα	19,21,24,25
α-MSH	22,25,26	β-2-Microglobulin	23,25,26,28	CD27	21,25
α-Synuclein	22	β-Endorphin	21,25,26	CD28	21
α-Synuclein (autoantibody)	22	β2-Glycoprotein (autoantibody)	19	CD30/TNFRSF8 (soluble)	19,21,24
AAV1 (serology)	20	β-Tubulin (total)	32	CD31/PECAM-1 (soluble)	21,22,25,27
AAV2 (serology)	20	<i>B. burgdorferi</i> Flagellin	20	CD34 (soluble)	27
AAV5 (serology)	20	<i>B. burgdorferi</i> OspA	20	CD40	21,25
AAV6 (serology)	20	<i>B. burgdorferi</i> OspB	20	CD40L (soluble)	17,18,21,23,24,25
AAV8 (serology)	20	<i>B. burgdorferi</i> OspC	20	CD44	22
AAV9 (serology)	20	<i>B. burgdorferi</i> VlsE1	20	CD63	21
AChR (autoantibody)	22	<i>B. microti</i> p41	20	CD80/B7-1	21,25
ACTH	17,21,23,24,26,27	B7-H2/ICOSL	21,25	CD81	21
Active Caspase 8 (total)	29	B7-H3/CD276	21,25	CD86/B7-2	21
Active Caspase 9 (total)	29	B7-H4/VTCN1	21	CD137/4-1BB/TNFRSF9	
ADAMTS13	21	B7-H5/VISTA	21	(soluble)	18,21,23,24,25
Adiponectin	20,24,25,26	B7-H6	21	CD226/DNAM-1	21,25
Adipsin/Factor D	19,20,21,23,25	BAD (total)	29	CD9	21
AGNA (SOX1) (autoantibody)	22	BAD (Ser112)	29	CEA	21
Agouti-Related Protein (AgRP)	21,23	BAFF (total)	18,19,21,23	CENP-A (autoantibody)	19
ALK-1 (soluble)	25	BAFF/Blys (autoantibody)	19	CENP-B (autoantibody)	19
Akt/PKB (Ser473)	29,30,31	BAX (total)	29	CHI3L1/YKL40	19,22,23
Akt/PKB (total)	29,30,31	BCA-1/CXCL13	17,18,23,25	CHIKV E	20
Akt1 (Ser473)	31	Bcl-2 (total)	29	Chk1 (Ser345)	30
Akt1 (total)	31	Bcl-2 (Ser70)	29	Chk2 (Thr68)	30
Akt2 (Ser474)	31	Bcl-xL (total)	29	CHOP (Total)	29
Akt2 (total)	31	Bcl-xL/BAD (interaction)	29	CK-MB	17,21
Albumin	23,26,28	BDNF	17,20,22,24,26	Cleaved PARP (total)	29
ALDH1A1	22	Betacellulin	23,25	Clusterin	23,25,26,28
AMPA -R (autoantibody)	22	BFV Polypeptide	20	CNTF	21,23
Amphiphysin (autoantibody)	22	BIM	29	Collagen IV	23
Amphiregulin	25	Blk (Tyr389)	31	COMP (cartilage oligomeric matrix protein)	19
Amylin (active)	17,20,23,24,26	BMP-9	22	Complement C1q	19,23
Amylin (total)	20,27	BNP	17,21	Complement C1ra	23
Amyloid beta 1-40 (Aβ 1-40)	22,25	BRAF (Ser446)	31	Complement C2	19
Amyloid beta 1-42 (Aβ 1-42)	22,25	BRAK/CXCL14	19	Complement C3	19,22
Amyloid beta 1-42 (Aβ 1-42) (autoantibody)	22	BTLA	21,25	Complement C3b/iC3b	19,23
<i>A. phagocytophilum</i> Msp5	20	C		Complement C4	19,22
<i>A. phagocytophilum</i> p44	20	c-Jun (Ser73)	30	Complement C4b	19,23
Angiogenin (ANG)	22	c-Kit (panTyr)	30	Complement C5	19
Angiopoietin-2	22,25	c-Kit/SCFR (soluble)	22	Complement C5a	19,23
Angiostatin/Kringle	22	c-Met/HGFR (panTyr)	30	Complement C5b-9	19
Angiotensinogen (AGT)	22	c-Myc (total)	30	Complement C9	19
ANGPTL3	20	C-Peptide	17,20,23	Complement Factor B	19
ANGPTL4	20	C-Peptide 2	24,26	Complement Factor H	19,22,23
ANGPTL6/AGF	20	C-Reactive Protein (CRP)	21,22,25,27	Complement Factor I	19
ANNA2 (RI)/NOVA1 (autoantibody)	22	C1q (autoantibody)	19	Complex IV (Total)	29
Apelin	20	CA125	21	Connective Tissue Growth Factor (CTGF)	26
Apolipoprotein AI (Apo AI)	21,22	CA15-3	21	Contactin-1	22
Apolipoprotein AII (Apo AII)	21	CA19-9	21	Corticosterone	26
Apolipoprotein B (Apo B)	21	Calbindin	23	Cortisol	22,28
Apolipoprotein CII (Apo CII)	21	Calreticulin (CALR)	21	CRAF (Ser338)	31
Apolipoprotein CIII (Apo CIII)	21,22	Carbonic Anhydrase 9 (CA9)	22	Creatine Kinase Muscle (CKM)	26
		Cardiac Troponin I (cTnI)	17,21,25,26	CREB (Ser133)	30,31

Analytes in blue font: available in Screening and Qualified assays. **Analytes in black font:** available in Qualified assays only.
Human pp. 17-23; Non-Human Primate p. 23; Mouse pp. 23-25; Rat pp. 25-26; Other/Multi-Species pp. 27-28; Cell Signaling kits pp. 29-32

MILLIPLEX® Analyte Index (continued)

Analyte	Page No.	Analyte	Page No.	Analyte	Page No.
CREB (total)	30,31	Ferritin	21,22	HCC-1/CCL14	19
CRMP5 (autoantibody)	22	Fetuin A	21,22	HCC-4/CCL16	19
CRYAB (autoantibody)	22	FGF-1/FGF-acidic	22	HCGβ	21
CTACK/CCL27	18,20,24	FGF-2/FGF-basic	17,18,21,22,23,25,27	Hck (Tyr411)	31
CTAG1B/NY-ESO-1 (autoantibody)	21	FGF-19	20	HE4	21
CTLA-4/CD152	21,25	FGF-21	17,20,24,26	Hepsin	22
CXCL6/GCP-2	17,18,19,21	FGF-23	17,20,21,26	HER2/ErbB2 (autoantibody)	21
CXCL16	18,19,21,23,25	FGFR1 (panTyr)	31	HGF	17,20,21,22,25
Cyclin B1 (Total)	29	FGL1/Hepassocin	21	HGFR/c-Met (soluble)	22
CYFRA21-1	21	Fgr (Tyr412)	31	HIF-1α (Total)	29
Cystatin C	23,25,26,28	Fibrinogen	21,27	HIF-1α (autoantibody)	21
D		Fibroblast Activation Protein (FAP)	21	Histone H3 (Ser10)	29
D-dimer	21	Flotillin-1 (FLOT-1)	21	HMGB1	18,19
DENV-1 NS1	20	Flt3 Ligand	17,18,23,27	HSP27 (Ser78)	30
DENV-2 NS1	20	Flt3 (panTyr)	30	HSP60 (autoantibody)	21
DENV-3 NS1	20	Follistatin (FST)	21,22,25	HSP70	19
DENV-4 NS1	20	Follistatin-like Protein 1 (FSTL1)	24,26	HSP70 (Total)	29
DKK1	17,21,22,26	Fractalkine/CX3CL1	17,18,20,23,25,26,27	Human Serum Albumin (HSA)	23
DPP4	20	<i>F. tularensis</i> LptD	20	HVEM	21,25
DPPX (autoantibody)	22	FSH	21,23,24,26,27	I	
DRD2 (autoantibody)	22	Fyn (Tyr420)	31	I-309/CCL1	17,18
E		G		ICAM-1 (soluble)	19,21,22,23,25,26
E-Cadherin	21,25	Galectin-1	21,25	ICOS	21
EEEV Spike E2	20	Galectin-1 (GAL1) (autoantibody)	21	IDO1	21
EGF	17,18,22,23,25,26,27	Galectin-3	21,25	IFNα	23,27
EGFR/HER1/ErbB1 (soluble)	19,22	Galectin-3 (GAL3) (autoantibody)	21	IFNα2	17,18,19,23
EGFR2/Her2/ErbB2 (soluble)	22	GAPDH (total)	21,32	IFNα2 (autoantibody)	19
EGFR3/Her3/ErbB3 (soluble)	22	G-CSF	17,18,22,23,25,27	IFNβ	18,19,23,27
ELAVL4 (HuD) (autoantibody)	22	G-CSF (autoantibody)	19	IFNβ (autoantibody)	19
eIF-4B (Ser422)	29,30	GDF-11	20	IFNγ	17,18,19,23,24,25,27,28
eIF-4E (Ser209)	30	GDF-15	20,21,22,23	IFNγ (autoantibody)	19
eIF-4G (Ser1108)	30	GDNF	22	IFNγ Receptor 1 (IFNGR1)	19
eIF2α (Ser51)	30	GH	21,23,24,26,27	IFNε	19
ENA-78/CXCL5	17,18	Ghrelin (active)	17,20,23,24,26,27	IFNω	18,19
Endocan-1 (ESM-1)	21,25	GIP (total)	17,20,23,24,26,27	IFNω (autoantibody)	19
Endoglin	22,25	GITR	21,25	IgA	20,24
Endostatin	21	GITRL	21	IgE	20,24
Endothelin-1	22,25	Glial Fibrillary Acidic Protein (GFAP) ...	22	IGF-1	20
ENO1 (autoantibody)	21	Glial Fibrillary Acidic Protein		IGF-2	20
Eotaxin/CCL11	17,18,23,25,27	(GFAP) (autoantibody)	22	IGF1R (total)	29
Eotaxin-2/CCL24/MP1F-2	17,18	GLP-1 (active)	20,23,24,26,27	IGF1R (panTyr)	30
Eotaxin-3/CCL26	17,18	GLP-1 (total)	17,20,24,26	IGF1R (Tyr1135/1136)	29
EpCAM	22	Glucagon	17,20,23,24,26,27	IGFBP1	20
Erythropoietin (EPO)	17,20,23,24,26	GM-CSF	17,18,23,24,25,27,28	IGFBP2	20
ErbB1/EGFR (panTyr)	30	GM-CSF (autoantibody)	19	IGF2BP2/IMP2 (autoantibody)	21
ErbB2/HER2 (panTyr)	30	GnRH	20,24	IGFBP3	20,21
ErbB3/HER3 (panTyr)	30	GOT1	26	IGF2BP3/IMP3 (autoantibody)	21
ErbB4/HER4 (panTyr)	30	GP130 (soluble)	19,24	IGFBP4	20
Erk/MAPK 1/2 (Thr185/Tyr187)	30,31	Granulysin	21	IGFBP5	20
Erk/MAPK 1/2 (total)	30,31	Granzyme A	17,18,23	IGFBP6	20
E-Selectin (soluble)	21,22,25,26	Granzyme B	17,18,19,21,23,24,25	IGFBP7	20
Estradiol	28	GROα/CXCL1	17,18,23,24,25,26,27	IgG1	20,24
Exodus-2/CCL21/6Ckine	17,18	GRP78 (Total)	29	IgG2	20
F		GSK3α (total)	29	IgG2a	24
FABP1	20,23	GSK3α (Ser21)	29	IgG2b	24
FABP3	17,20,21,22,26	GSK3β (total)	29	IgG3	20,24
FABP4	17,21	GSK3β (Ser9)	29	IgG4	20
FADD (Ser194)	30	GSTα	23,26	IgM	20,24
Fas Ligand (sFasL)		H		IκBα (Ser32)	30
(soluble)	17,18,19,21,23,24,27	H2AX (Ser139)	29,30	IKKα/β (Ser177/181)	30
Fas/TNFRSF6		Haptoglobin	21,25,26,27	IL-1α	17,18,19,23,24,25,27
(soluble)	17,18,19,21,23,24,25	HB-EGF	17,22	IL-1α (autoantibody)	19
				IL-1β	17,18,20,21,23,24,25,27,28
				IL-1RI (soluble)	19,24

Analytes in blue font: available in Screening and Qualified assays. **Analytes in black font:** available in Qualified assays only.

Human pp. 17-23; Non-Human Primate p. 23; Mouse pp. 23-25; Rat pp. 25-26; Other/Multi-Species pp. 27-28; Cell Signaling kits pp. 29-32

MILLIPLEX® Analyte Index (continued)

Analyte	Page No.	Analyte	Page No.	Analyte	Page No.
IL-1RII (soluble).....	19,24	IRS1 (Ser636).....	29	MICA.....	21
IL-1Ra.....	17,18,23,27,28	IRS1 (total).....	29	MICB.....	21
IL-2.....	17,18,23,24,25,27,28	I-TAC/CXCL11.....	17,18,19,23	Midkine.....	22
IL-2 (autoantibody).....	19	J		MIF.....	19,21
IL-2Ra (soluble).....	19,21,24,25	Jag1.....	20,24	MIG/CXCL9.....	17,18,19,23,25,27
IL-3.....	17,18,23,25	JEV E.....	20	MIP-1a/CCL3.....	17,18,19,23,25,27
IL-3 (autoantibody).....	19	JEV NS1.....	20	MIP-1β/CCL4.....	17,18,19,23,24,27
IL-4.....	17,18,23,24,25,27,28	JNK/SAPK1 (Thr183/Tyr185).....	29,30,32	MIP-1δ/MIP-5/CCL15.....	17,18
IL-4 (autoantibody).....	19	JNK/SAPK1 (total).....	30,32	MIP-2/CXCL2.....	23,24,25
IL-4R (soluble).....	19,24	Jo-1 (autoantibody).....	19	MIP-3a/CCL20.....	17,18,19,23,24,27
IL-5.....	17,18,23,24,25,27	K		MIP-3β/CCL19.....	17,18,19,23
IL-6.....	17,18,20,21,23,24,25,26,27,28	Kallikrein-6.....	22	MIP-4/PARC/CCL18.....	19,22
IL-6 (autoantibody).....	19	KC-like Protein.....	27	MMP-1.....	17,19,27
IL-6R (soluble).....	19,24	KC/GRO/CINC-1/CXCL1.....	23,24,25,26,27	MMP-2.....	17,19,24
IL-6Ra (soluble).....	22	KIM-1.....	23,25,26,28	MMP-3.....	17,19,24
IL-7.....	17,18,23,24,25,27	Ki-67 (Total).....	29	MMP-7.....	17,19
IL-8/CXCL8.....	17,18,19,20,21,22,23,25,27,28	Ku (autoantibody).....	19	MMP-8.....	19,24
IL-8/CXCL8 (autoantibody).....	19	L		MMP-9.....	17,19,27
IL-9.....	17,18,23	Lactotransferrin (LTF).....	19	MMP-10.....	17,19
IL-10.....	17,18,20,23,24,25,27	LAG-3.....	21,25	MMP-12.....	17,19,24
IL-10 (autoantibody).....	19	LAT (panTyr).....	31	MMP-13.....	19
IL-11.....	17,18,19,23	LC3B (Total).....	29	MOG (autoantibody).....	22
IL-12.....	27,28	Lck (panTyr).....	31	MPIF/CCL23.....	18,19
IL-12 (p40).....	17,18,23,25,27	Lck (Tyr394).....	31	MSK1 (Ser212).....	30
IL-12 (p40) (autoantibody).....	19	Leishmania K39.....	20	mTOR (Ser2448).....	29,32
IL-12 (p70).....	17,18,23,24,25,27	Leptin.....	17,20,21,22,23,24,25,26,27	mTOR (total).....	29,32
IL-13.....	17,18,23,24,25,27	Leucine-rich alpha-2-glycoprotein 1 (LRG1).....	21	Mucin-1 (MUC1) (autoantibody).....	21
IL-14/α -Taxilin.....	19	Leukemia Inhibitory Factor (LIF).....	17,18,20,23,24,26	Myeloperoxidase (MPO).....	21,22
IL-15.....	17,18,20,23,24,25,26,27	LH.....	21,23,24,26,27	Myeloperoxidase (MPO) (autoantibody).....	19
IL-15 (autoantibody).....	19	LIGHT.....	21,25	Myoglobin.....	21
IL-16.....	17,18,23,24,25,27	LIX/CXCL5.....	23,24,25	Myosin Light Chain 3 (MYL3).....	26
IL-17A/CTLA8.....	17,18,23,24,25,27	L-Selectin (soluble).....	21	Myostatin/GDF8.....	17,20,24,26
IL-17A/CTLA8 (autoantibody).....	19	Lymphotactin/XCL1.....	18,19	N	
IL-17E/IL-25.....	17,18,19,23,24	Lyn (Tyr397).....	31	NAP-2/CXCL7.....	19
IL-17F.....	17,18,23,24	M		NCAM (soluble).....	22
IL-17F (autoantibody).....	19	Mannose-Binding Lectin (MBL).....	19,23	NCAM1/L1CAM/CD171.....	22
IL-18.....	17,18,19,20,23,25,27,28	MAYV E.....	20	Nectin-2.....	21
IL-18 (autoantibody).....	19	MBP (autoantibody).....	22	Nectin-4.....	21
IL-18BPα.....	19	MCAM/MUC18/sCD146.....	21	Neurogranin (NRGN).....	22
IL-19.....	19	Mcl-1/BIM (interaction).....	29	Neuron-specific Enolase (NSE).....	22
IL-20.....	18,23	Mcl-1 (total).....	29	Neuropilin-1 (soluble).....	22
IL-21.....	17,18,23,24,27	MCP-1/CCL2.....	17,18,20,23,24,25,26,27,28	Neurotensin.....	22,25,26
IL-22.....	17,18,23,24	MCP-2/CCL8.....	17,18,23	Neutrophil Elastase-2 (ELA2).....	19
IL-22 (autoantibody).....	19	MCP-3/CCL7.....	17,18,23	NFκB (total).....	30
IL-23.....	17,18,23,24	MCP-4/CCL13.....	17,18	NFκB (Ser536).....	30
IL-24.....	18,19	MCP-5/CCL12.....	23	NF-L (autoantibody).....	22
IL-27.....	18,24	M-CSF.....	17,18,19,23,25,27	Nipah Virus gpG.....	20
IL-28A/IFNλ2.....	18,19,23	M-CSFR (panTyr).....	30	NGAL/Lipocalin-2.....	19,20,21,23,25,26,28
IL-28B/IFNλ3.....	19,23,24,27	MDC/CCL22.....	17,18,23	NGF.....	17,20
IL-29/IFNλ1.....	18,19	MDH1.....	23	NMDAR1 (autoantibody).....	22
IL-31.....	18,23,24	MDM2 (total).....	30	Notch1.....	24
IL-32α.....	19	MEK1 (Ser217/221).....	31	NOXA/Mcl-1 (interaction).....	29
IL-33/NF-HEV (mature).....	17,18,23,24	MEK1 (Ser222).....	30	NRXN3 (autoantibody).....	22
IL-34.....	18,19	Melanoma Inhibitory Activity Protein (MIA).....	21	NT-proBNP.....	17,21
IL-35.....	18,19	Melatonin.....	22,26	O	
IL-36β/IL-1F8.....	19	Mesothelin.....	22	Oncostatin-M (OSM).....	20,21,24,25
IL-36Ra/IL-1F5.....	27	mGluR1 (autoantibody).....	22	OROV Nucleoprotein.....	20
IL-37/IL-1F7.....	19	Mi-2 (autoantibody).....	19	Osteoactivin.....	23
IL-38/IL-1F10.....	19			Osteocalcin (OC).....	17,21
Insulin.....	17,20,21,23,24,25,26,27			Osteocarin (OSTN)/Musclin.....	20,24,26
IP-10/CXCL10.....	17,18,23,25,26,27			Osteonectin/SPARC.....	17,20,22,24,26
IR (total).....	29			Osteopontin (OPN).....	17,21,22,23,25,26,28
IR (Tyr162/1163).....	29				
Irisin.....	20,24,26				

Analytes in blue font: available in Screening and Qualified assays. **Analytes in black font:** available in Qualified assays only.

Human pp. 17-23; Non-Human Primate p. 23; Mouse pp. 23-25; Rat pp. 25-26; Other/Multi-Species pp. 27-28; Cell Signaling kits pp. 29-32

MILLIPLEX® Analyte Index (continued)

Analyte	Page No.	Analyte	Page No.	Analyte	Page No.
Osteopontin (OPN) (autoantibody)	19	Ras (total).....	31	Tau (Total)	22
Osteoprotegerin (OPG)	17,21,22,26	RasG12D (total)	31	Tau (Total) (autoantibody)	22
OX40/CD134.....	21	RasG12V (total)	31	Tau (Thr181)	22
Oxytocin	22,25,26	RBP4	23,28	TBEV NS1 (Europe)	20
P		Renin	23,25	TBEV NS1 (Far East).....	20
P16-INK4a (autoantibody)	21	Resistin	19,20,24	Testosterone	28
p21 (total).....	30	Ribosomal P (autoantibody).....	19	TFF-3.....	23,28
p38/SAPK2A/B (Thr180/Tyr182)	30,32	<i>R. rickettsii</i> Omp.....	20	TGFa	17,18,21,23,27
p38/SAPK2A/B (total).....	30,32	<i>R. typhi</i> Omp.....	20	TGFβ1.....	28
p53 (autoantibody).....	21	RNA Pol II (Ser 2)	29	TGFβ2.....	28
p53 (Ser15).....	29,30	RNP (autoantibody)	19	TGFβ3.....	28
p53 (Ser46).....	29	RNP/Sm (autoantibody)	19	TGFβRII (total)	31
p70S6 Kinase (total).....	29,30	RPS6 (Ser235/Ser236).....	29	Thrombomodulin.....	21,25
p70S6 Kinase (Thr389/412).....	29,30	RPS6 (total)	29	Thrombospondin-1 (TSP-1)	19
PAI-1 (total)	19,20,22,24,25,26,27	RRV Polyprotein	20	Thrombospondin-2 (TSP-2)	22
Pancreatic Polypeptide (PP)		RVFV Nucleoprotein	20	TIE1 (panTyr)	30
.....	17,20,23,24,26,27	S		TIE2 (panTyr)	30
PARK5/UCHL1.....	22	S100B.....	22	Tie-2 (soluble)	22
PARK7/DJ1	22	SARS-CoV-2 N	20	TIM-3	21,25
PCA1 (Yo)/CDR2 (autoantibody).....	22	SARS-CoV-2 RBD	20	TIMP-1	17,19,23,25,26
PCNA (autoantibody)	19	SARS-CoV-2 Spike S1	20	TIMP-2.....	19
PD-1	21,25	SARS-CoV-2 Spike S2	20	TIMP-3.....	19
PDGF-AA	17,18,22	SCF	17,18,21,27	TIMP-4	19
PDGF-AB/BB	17,18,22	Scl-70 (autoantibody)	19	Tissue Factor (TF).....	21
PDGF-BB.....	27	Sclerostin (SOST)	17,21,26	TLR-2	21,25
PDGFRβ (panTyr)	30	SDF-1/CXCL12	17,18,25,27	TNFR1 (total)	30
PD-L1	21,25	SDH	23,26	TNFa	17,18,20,21,23,24,25,26,27,28
PD-L2	21,25	Secretin.....	20,24,26	TNFa (autoantibody).....	19
PEDF	22	Serum Amyloid A (SAA)	21	TNFB/Lymphotoxin-α (LTA)	17,18,23,24
Pentraxin-3 (PTX3).....	21	Serum Amyloid P (SAP).....	21,22,25	TNF RI (soluble) (sTNFR1)	19,24
Perforin	17,18,21,23	Sex Hormone Binding Globulin (SHBG) 21		TNF RII (soluble) (sTNFR2)	19,24
Periostin	22	Siglec-7	21	TPO	17,18
Polyubiquitin K48-linkage (Total)	29	Siglec-9	21	TRAIL/TNFSF10	17,18,21
Progranulin (pGRN)	21	SLEV NS1.....	20	Transglutaminase 2 (TGM2)	22
PL-12 (Alanyl-tRNA Synthetase)		Sm (autoantibody)	19	TRAP5	22
(autoantibody).....	19	SMAD2 (Ser465/467).....	31	TREM2	22
Placental Growth Factor 2	25	SMAD3 (Ser423/425).....	31	Troponin I (TnI)	17,21,25,26
Placental Growth Factor (PIGF)		SMAD4 (total)	31	Troponin T (TnT)	21,25,26
.....	17,21,22,25	SOX2 (autoantibody)	21	<i>T. cruzi</i> /Chagas.....	20
<i>P. falciparum</i> LDH	20	Src (Tyr419)	31	TSC2 (total)	29
<i>P. vivax</i> LDH.....	20	SSA/Ro52 (autoantibody).....	19	TSC2 (Ser939).....	29
Platelet Factor 4 (PF4)/CXCL4.....	21	SSA/Ro60 (autoantibody)	19	TSG101.....	21
PF4 (autoantibody).....	19	SSB/La (autoantibody)	19	TSH.....	21,23,24,26,27
PM/Sci 100 (autoantibody)	19	STAT1 (Tyr701)	30,31	TSLP	17,18
PNMA1 (Ma1) (autoantibody)	22	STAT2 (Tyr690)	31	TWEAK	22
POWV NS1	20	STAT3 (total).....	30,32	U	
Progesterone.....	28	STAT3 (Ser727).....	30	ucMGP	21
Prolactin	21,24,25,26,27	STAT3 (Tyr705).....	32	uPAR (soluble)/CD87	22
ProMMP-9	24,25	STAT5A/B (total).....	30	Uromodulin	23
Properdin	19,23	STAT5A/B (Tyr694/699).....	30,31	Uteroglobin	19
Proteinase 3 (autoantibody)	19	STAT6 (Tyr641).....	31	V	
PSA (free)	21	Substance P	22,25,26	VCAM-1 (soluble)	19,21,22
PSA (total)	21	Superoxide Dismutase 1 (soluble)		VEGF	25,26,27
P-Selectin (soluble)	21,25	(sSOD1)	22	VEGF-A	17,18,21,22,23,25,26,27
PTEN (total)	29	Superoxide Dismutase 2 (soluble)		VEGF-C	17,22,25
PTEN (Ser380)	29	(sSOD2)	22	VEGF-D	17,22,25
PTH	21,23,26	Survivin/BIRC5 (autoantibody)	21	VEGFR1/Flt-1 (panTyr)	30
PVR/CD155	21	Syntenin-1 (Syn-1)	21	VEGFR1/Flt-1 (soluble).....	19,22,24
PYY (total)	17,20,23,24,26,27	Syk (panTyr).....	31	VEGFR2/KDR/Flk-1 (panTyr)	30
R		T		VEGFR2/KDR/Flk-1 (soluble).....	19,22,24
RAGE (soluble)	19,22,23,24	T3	26,28	VEGFR3/Flt-4 (panTyr).....	30
RANKL	21	T4	26,28	VEGFR3/Flt-4 (soluble)	19,22,24
RANTES/CCL5	17,18,22,23,24,25,27	TARC/CCL17	17,18,23	Vitronectin	23

Analytes in blue font: available in Screening and Qualified assays. **Analytes in black font:** available in Qualified assays only.
Human pp. 17-23; Non-Human Primate p. 23; Mouse pp. 23-25; Rat pp. 25-26; Other/Multi-Species pp. 27-28; Cell Signaling kits pp. 29-32

W

WEEV Spike E1	20
WNV E	20
WNV NS1	20

Y

YFV E	20
YFV NS1	20
Yes (Tyr421)	31
YKL40/CHI3L1	19,22,23

Z

ZAP-70 (panTyr)	31
ZIKV E	20
ZIKV NS1 (French Polynesia)	20
ZIKV NS1 (Uganda)	20

Analytes in blue font: available in Screening and Qualified assays. **Analytes in black font:** available in Qualified assays only.

Human pp. 17-23; Non-Human Primate p. 23; Mouse pp. 23-25; Rat pp. 25-26; Other/Multi-Species pp. 27-28; Cell Signaling kits pp. 29-32

What's In A Typical Qualified MILLIPLEX® Kit?



Multiplex Your Research

Explore all the ways researchers are using MILLIPLEX® multiplex assays at SigmaAldrich.com/milliplex-learn



MILLIPLEX® PLEXpedition Screening Assays

HUMAN

Multiple Research Areas

Human Screening Assay

(Cat. No. HPLX1-115SP)

80 (Cat. No. HPLX1-115SP-PX80) ♦

115 (Cat. No. HPLX1-115SP-PX)

xMAP® INTELLIFLEX® or FLEXMAP 3D® Instrument						
Luminex® 200™ Instrument						
MAGPIX® Instrument						
EGF♦	IL-2♦	MCP-1♦	ACTH♦	M-CSF♦	Amylin (Active)	MCP-2
FGF-2♦	IL-4♦	MIG♦	BDNF♦	MCP-3♦	BCA-1	MCP-4
Fractalkine♦	IL-5♦	MIP-1α♦	C-Peptide♦	MDC♦	BNP	MIP-16
G-CSF♦	IL-6♦	MIP-1β♦	Eotaxin♦	MIP-3α♦	CK-MB	MIP-3β
Ghrelin (Active)♦	IL-7♦	OPN♦	FABP4♦	MMP-1♦	GCP-2	MMP-3
GIP♦	IL-8♦	PDGF-AA♦	FLT3 Ligand♦	MMP-2♦	DKK1	MMP-10
Glucagon♦	IL-10♦	PDGF-AB/BB♦	GLP-1 Total♦	MMP-7♦	ENA-78	MMP-12
GM-CSF♦	IL-12 (p40)♦	Perforin♦	Granzyme A♦	MMP-9♦	Eotaxin-2	Myostatin
Granzyme B♦	IL-12 (p70)♦	RANTES♦	I-TAC♦	OPG♦	Eotaxin-3	NGF
GROα♦	IL-13♦	sCD40L♦	IL-3♦	PIGF♦	Erythropoietin (EPO)	NTproBNP
HGF♦	IL-15♦	SCF♦	IL-9♦	PP♦	FABP3	Osteocalcin
IFNa2♦	IL-17A♦	sFAS♦	IL-17F♦	PYY♦	FGF-21	Osteonectin
IFNγ♦	IL-18♦	sFASL♦	IL-21♦	TARC♦	FGF-23	SDF-1
IL-1α♦	IL-22♦	TGFα♦	IL-33♦	TNFβ♦	FGF-23	SOST
IL-1β♦	IL-23♦	TNFα♦	Insulin♦	Troponin I♦	HB-EGF	TPO
IL-1Ra♦	IP-10♦	TRAIL♦			I-309	TSLP
	Leptin♦	VEGF-A♦			IL-11	VEGF-C
					LIF	VEGF-D

Configure Your Multiplex Kits with Ease

Check out our enhanced kit configurator where you can select your analytes for MILLIPLEX® qualified assays as well as our new MILLIPLEX® PLEXpedition screening kits. Design your kit at SigmaAldrich.com/milliplex-order



MILLIPLEX® Qualified Assays

HUMAN

Immunology

Human Cytokine/Chemokine/ Growth Factor Panel A

☐ (Cat. No. HCYTA-60K)

38 (Cat. No. HCYTA-60K-PX38)♦

38 (Cat. No. HCYTA-60K-PXBK38)♦

48 (Cat. No. HCYTA-60K-PX48)

48 (Cat. No. HCYTA-60K-PXBK48)

sCD40L	IL-13♦
EGF♦	IL-15♦
Eotaxin/CCL11♦	IL-17A/CTLA8♦
FGF-2/FGF-basic	IL-17E/IL-25♦
Flt3 Ligand	IL-17F♦
Fractalkine/CX3CL1	IL-18♦
G-CSF♦	IL-22♦
GM-CSF♦	IL-27
GROα	IP-10/CXCL10♦
IFNα2♦	MCP-1/CCL2♦
IFNγ♦	MCP-3/CCL7
IL-1α♦	M-CSF♦
IL-1β♦	MDC/CCL22
IL-1RA♦	MIG/CXCL9♦
IL-2♦	MIP-1α/CCL3♦
IL-3♦	MIP-1β/CCL4♦
IL-4♦	PDGF-AA♦
IL-5♦	PDGF-AB/BB♦
IL-6♦	RANTES/CCL5♦▲*
IL-7♦	TGFα
IL-8/CXCL8♦	TNFα♦
IL-9	TNFβ/Lymphotoxin-α (LTA)♦
IL-10♦	VEGF-A♦
IL-12 (p40)♦	
IL-12 (p70)♦	

Get both Human Panel A and Panel B premixed kits with one catalog number!

76-Plex (Two 38-plex premixed kits)

Cat. No. HCYTPAB-76SK

96-Plex (Two 48-plex premixed kits)

Cat. No. HCYTPAB-96SK

Kits cannot be combined.

Human Cytokine/Chemokine/ Growth Factor Panel B

☐ (Cat. No. HCYTB-60K)

38 (Cat. No. HCYTB-60K-PX38)♦

38 (Cat. No. HCYTB-60K-PXBK38)♦

48 (Cat. No. HCYTB-60K-PX48)

48 (Cat. No. HCYTB-60K-PXBK48)

6Ckine/CCL21♦	IL-28A/IFNA2♦
APRIL♦	IL-29/IFNA1♦
BAFF/Blys♦	IL-31♦
BCA-1/CXCL13♦	IL-33/NF-HEV (mature)♦
CCL28	IL-34
sCD137/4-1BB/TNFRSF9♦	IL-35
CTACK/CCL27♦	IFNβ♦
CXCL6/GCP-2	IFNω♦
CXCL16	I-TAC/CXCL11♦
ENA-78/CXCL5♦	LIF♦
Eotaxin-2/CCL24/MIPF-2♦	Lymphotoxin/XCL1
Eotaxin-3/CCL26♦	MCP-2/CCL8♦
sFAS/TNFRSF6♦	MCP-4/CCL13♦
sFASL♦	MIP-1δ/MIP-5/CCL15♦
Granzyme A♦	MIP-3α/CCL20♦
Granzyme B♦	MIP-3β/CCL19
HMGB1♦	MIPF-1/CCL23
I-309/CCL1♦	Perforin♦
IL-11	SCF♦
IL-16♦	SDF-1/CXCL12♦
IL-20♦	TARC/CCL17♦
IL-21♦	TPO♦
IL-23♦	TRAIL/TNFSF10♦
IL-24	TSLP♦

Human High Sensitivity T Cell

☐ (Cat. No. HSTCMAG-28SK)

13 (Cat. No. HSTCMAG28SPMX13)♦

13 (Bulk Cat. No. HSTCMAG28SPMX13BK)♦

21 (Cat. No. HSTCMAG28SPMX21)

21 (Bulk Cat. No. HSTCMAG28SPMX21BK)

Fractalkine/CX3CL1	IL-12 (p70)♦
GM-CSF♦	IL-13♦
IFNγ♦	IL-17A/CTLA8
IL-1β♦	IL-21
IL-2♦	IL-23
IL-4♦	I-TAC/CXCL11
IL-5♦	MIP-1α/CCL3
IL-6♦	MIP-1β/CCL4
IL-7♦	MIP-3α/CCL20
IL-8/CXCL8♦	TNFα♦
IL-10♦	

384-Well Human High Sensitivity T Cell

(Cat. No. HSTC384-28K)

21 (Cat. No. HSTCMAG384-PX21)

21 (Bulk Cat. No. HSTCMAG384PX21BK)

Fractalkine/CX3CL1	IL-12 (p70)
GM-CSF	IL-13
IFNγ	IL-17A/CTLA8
IL-1β	IL-21
IL-2	IL-23
IL-4	I-TAC/CXCL11
IL-5	MIP-1α/CCL3
IL-6	MIP-1β/CCL4
IL-7	MIP-3α/CCL20
IL-8/CXCL8	TNFα
IL-10	

Human Th17

(Cat. No. HTH17MAG-14K)

25 (Cat. No. HT17MG-14K-PX25)

25 (Bulk Cat. No. HT17MAG14PMX25BK)

GM-CSF	IL-17E/IL-25
IFNγ	IL-17F
IL-1β	IL-21
IL-2	IL-22
IL-4	IL-23
IL-5	IL-27
IL-6	IL-28A/IFNA2
IL-9	IL-31
IL-10	IL-33/NF-HEV (mature)
IL-12 (p70)	MIP-3α/CCL20
IL-13	TNFα
IL-15	TNFβ/Lymphotoxin-α (LTA)
IL-17A/CTLA8	

Human CD8+ T Cell

☐ (Cat. No. HCD8MAG-15K)

17 (Cat. No. HCD8MAG15K17PMX)

17 (Bulk Cat. No. HCD8MAG15KPX17BK)

sCD137/4-1BB/TNFRSF9	IL-4
sFAS/TNFRSF6	IL-5
sFasL	IL-6
GM-CSF	IL-10
Granzyme A	IL-13
Granzyme B	MIP-1α/CCL3
IFNγ	MIP-1β/CCL4
IL-2	Perforin
	TNFα

MILLIPLEX® Immunology Kits

Profile more immune biomarkers faster with MILLIPLEX® multiplex immunoassays. Learn more at SigmaAldrich.com/milliplex-immunology



Legend key for MILLIPLEX® kits

† Analytes which cannot be plexed together:
– Active and total
– Free and total

♦ Available in Cat. No. listed

▼ Premix panel only

Bulk indicates Space Saver Packaging

▲ These analytes cannot be plexed with other analytes in this panel in serum/plasma

▷ Requires a protease inhibitor during sample collection

■ Requires sample extraction

** No standard/QC; antibody levels are reported as MFIs (Median Fluorescence Intensities)

● Serum/Plasma only

★ Tissue Culture samples only

^ Competitive assay format

☐ Available for custom premix

◦ Analyte requires a different sample dilution from others in panel

⊙ Analyte cannot be measured in serum

Requires sample acidification

§ Analyte cannot be measured in plasma

u Urine only

Human Interferon

(Cat. No. HIFN-130K)

9 (Cat. No. HIFN-130K-PMX9)

IFNα2	IFNε
IFNβ	IFNω
IFNγ	IL-28A/IFNλ2
IFNγ Receptor 1 (IFNGR1)	IL-28B/IFNλ3
	IL-29/IFNλ1

Human Soluble Cytokine Receptor

(Cat. No. HSCRMAG-32K)

14 (Cat. No. HSCRMAG32KPX14)

14 (Bulk Cat. No. HSCRMAG32PMX14BK)

sCD30	sIL-6R
sEGFR/sHER1/ErbB1	sRAGE
sGP130	sTNF RI
sIL-1RI	sTNF RII
sIL-1RII	sVEGFR1/sFlt-1
sIL-2Ra	sVEGFR2/sKDR/sFlk-1
sIL-4R	sVEGFR3/sFlt-4

Human Cytokine/Chemokine Panel III

(Cat. No. HCYP3MAG-63K)

HCC-1/CCL14▲	M-CSF
IL-11	MIG/CXCL9
IL-29/IFNλ1	MIP-3α/CCL20
I-TAC/CXCL11	MIP-3β/CCL19
CXCL6/GCP-2	NAP-2/CXCL7▲
Lymphotactin/XCL1	

Human Cytokine/Chemokine Panel IV

20 (Cat. No. HCYP4MAG-64K)

20 (Cat. No. HCY4MG-64K-PX20)

20 (Bulk Cat. No. HCP4MG64KPX20BK)

BAFF/Blys	IL-28B/IFNλ3
BRACK/CXCL14	IL-32α
CCL28	IL-34
CXCL16	IL-35
HCC-4/CCL16	IL-36β/IL-1F8
HMGB1●	IL-37/IL-1F7
IFNβ	IL-38/IL-1F10
IL-14/α-Taxilin	MIP-4/PARC/CCL18
IL-19	MPIF/CCL23
IL-24	YKL40/CHI3L1

Human IL-18 Singleplex

(Cat. No. HIL18MAG-66K)

IL-18

Human Complement Expanded Panel 1

10 (Cat. No. HCMPEX1-19K)

Adipsin/Factor D	C9
C2	Factor I
C4b	Mannose-Binding Lectin (MBL)
C5	Properdin
C5a	

Human Complement Panel 2

10 (Cat. No. HCMP2MAG-19K)

C1q	C4
C3	Factor B
C3b	Factor H

Human Complement C5b-9 Singleplex

(Cat. No. HC5B9-140K)

C5b-9

Human MMP Panel 1

(Cat. No. HMMPM1MAG-55K)

MMP-3	MMP-13
MMP-12	

Human MMP Panel 2

(Cat. No. HMMPM2MAG-55K)

MMP-1	MMP-9
MMP-2	MMP-10
MMP-7	

Human TIMP Panel 1 (Serum/Plasma samples)

(Cat. No. HTMP1MAG-54K)

TIMP-1	TIMP-2
--------	--------

Human TIMP Panel 2 (Cell culture samples)

(Cat. No. HTMP2MAG-54K)

TIMP-1	TIMP-3
TIMP-2	TIMP-4

Human Fibrosis Panel 1

(Cat. No. HFIB1-100K)

COMP	Uteroglobin
IL-18BPα	

Immune Response

Human Autoimmune Autoantibody IgG**

(Cat. No. HAIAB-10K)

QUALITATIVE KIT

β2-Glycoprotein	PM/ScI 100
C1q	Proteinase 3
CENP-A	Ribosomal P
CENP-B	RNP
Jo-1	RNP/Sm
Ku	ScI-70
Mi-2	Sm
Myeloperoxidase (MPO)	SSA/Ro52
PCNA	SSA/Ro60
PL-12 (Alanyl-tRNA Synthetase)	SSB/La

Human Cytokine Autoantibody Expanded IgG, IgA, and IgM **

IgG (Cat. No. HCYTABG-17K)

IgA (Cat. No. HCYTABA-17K)

IgM (Cat. No. HCYTABM-17K)

QUALITATIVE KIT

BAFF/Blys	IL-8/CXCL8
G-CSF	IL-10
GM-CSF	IL-12 (p40)
IFNα2	IL-15
IFNβ	IL-17A/CTLA8
IFNγ	IL-17F
IFNω	IL-18
IL-1α	IL-22
IL-2	Osteopontin
IL-3	PF4
IL-4	TNFα
IL-6	

Human Sepsis Panel 1

(Cat. No. HSP1MAG-63K)

sFAS/TNFRSF6	MIF
sFasL	PAI-1 (total)
sICAM-1	sVCAM-1

Human Sepsis Panel 2

(Cat. No. HSP2MAG-63K)

Granzyme B	MIP-1α/CCL3
HSP70	MIP-1β/CCL4
IL-1α	MMP-8
IL-8/CXCL8	

Human Sepsis Panel 3

(Cat. No. HSP3MAG-63K)

Lactotransferrin (LTF)	Resistin
Neutrophil Elastase-2 (ELA2)	Thrombospondin-1 (TSP-1)
NGAL/Lipocalin-2	

MILLIPLEX® Autoantibody Profiling

Discover the benefits of multiplexing autoantibodies at
SigmaAldrich.com/milliplex-autoab



Human Immunoglobulin Isotyping

(Cat. No. HGAMMAG-301K)

IgA	IgG3
IgG1	IgG4
IgG2	IgM

Human IgE – Singleplex

(Cat. No. HGAMMAG-303E)

IgE

Virology**Vector-Borne Disease Panel IgG**** **NEW**

(Cat. No. HVBBDG-270K)

- 6 Alphavirus (Cat. No. HVBBDG-270K-AVPMX) ◇
- 17 Flavivirus (Cat. No. HVBBDG-270K-FVPMX) ◆
- 42 (Cat. No. HVBBDG-270K-PMX)

QUALITATIVE KIT

<i>A. phagocytophilum</i> Msp5	MAYV E ◇
<i>A. phagocytophilum</i> p44	Nipah Virus gpG
<i>B. microti</i> p41	OROV Nucleoprotein
BFV Polyprotein ◇	<i>P. falciparum</i> LDH
<i>B. burgdorferi</i> Flagellin	<i>P. vivax</i> LDH
<i>B. burgdorferi</i> OspA	POWV NS1 ◆
<i>B. burgdorferi</i> OspB	<i>R. rickettsii</i> Omp
<i>B. burgdorferi</i> OspC	<i>R. typhi</i> Omp
<i>B. burgdorferi</i> VlsE1	RVFV Nucleoprotein
CHIKV E ◇	RRV Polyprotein ◇
CCHFV gp38	SLEV NS1 ◆
DENV-1 NS1 ◆	TBEV NS1 (Europe) ◆
DENV-2 NS1 ◆	TBEV NS1 (Far East) ◆
DENV-3 NS1 ◆	<i>T. cruzi</i> /Chagas
DENV-4 NS1 ◆	WNV E ◆
EEEV Spike E2 ◇	WNV NS1 ◆
<i>F. tularensis</i> LptD	WEEV Spike E1 ◇
JEV E ◆	YFV E ◆
JEV NS1 ◆	YFV NS1 ◆
Leishmania K39	ZIKV E ◆
	ZIKV NS1 (Uganda) ◆
	ZIKV NS1 (French Polynesia) ◆

AAV Serotyping Panel IgG**

(Cat. No. AAVSERG-230K)

QUALITATIVE KIT

AAV1	AAV6
AAV2	AAV8
AAV5	AAV9

SARS-CoV-2 Antigen Panel 1 IgA**

(Cat. No. HC19SERA1-85K)

QUALITATIVE KIT

SARS-CoV-2 N	SARS-CoV-2 Spike S1
SARS-CoV-2 RBD	SARS-CoV-2 Spike S2

SARS-CoV-2 Antigen Panel 1 IgG**

(Cat. No. HC19SERG1-85K)

QUALITATIVE KIT

SARS-CoV-2 N	SARS-CoV-2 Spike S1
SARS-CoV-2 RBD	SARS-CoV-2 Spike S2

SARS-CoV-2 Antigen Panel 1 IgM**

(Cat. No. HC19SERM1-85K)

QUALITATIVE KIT

SARS-CoV-2 N	SARS-CoV-2 Spike S1
SARS-CoV-2 RBD	SARS-CoV-2 Spike S2

For Research Use Only. Not For Use In Diagnostic Procedures.**Metabolism/Endocrinology****Human Adipokine Panel 1 (Serum/Plasma samples)**

(Cat. No. HADK1MAG-61K)

Adiponectin	PAI-1 (total)
Adipsin/Factor D	Resistin
NGAL/Lipocalin-2	

Human Adipokine Panel 2 (Serum/Plasma samples)

(Cat. No. HADK2MAG-61K)

IL-1β	Leptin
IL-6	MCP-1/CCL2
IL-8/CXCL8	NGF
Insulin	TNFα

Human Adipocyte (Cell culture samples)

(Cat. No. HADCYMAG-61K)

Adiponectin	MCP-1/CCL2
HGF	NGF
IL-1β	PAI-1 (total)
IL-6	Resistin
IL-8/CXCL8	TNFα
Leptin	

Human Metabolic Hormone V3

(Cat. No. HMH3-34K)

Amylin (active) †	Insulin
Amylin (total) †	Leptin
C-Peptide	MCP-1/CCL2
Ghrelin (active) †	Pancreatic Polypeptide (PP)
GIP (total)	PYY (total)
GLP-1 (active) †	Secretin
GLP-1 (total) †	TNFα
Glucagon †	
IL-6	

Human Myokine

(Cat. No. HMYOMAG-56K)

Apelin	IL-15
BDNF	Irisin
Erythropoietin (EPO)	LIF
FABP3	Myostatin/GDF8
FGF-21	Oncostatin-M (OSM)
Fractalkine/CX3CL1	Osteocrin (OSTN)/Musclin
FSTL1	Osteonectin/SPARC
IL-6	

Human Aging Panel 1

(Cat. No. HAGE1MAG-20K)

CTACK/CCL27	IL-6
FGF-21	IL-10*
GDF-11	IL-18
GDF15	Jag1
GnRH	Leptin

Human Liver Protein

(Cat. No. HLPPMAG-57K)

α-Fetoprotein (AFP)	FGF-19
ANGPTL3	FGF-21
ANGPTL4	FGF-23
ANGPTL6	HGF
FABP1	

Human IGF Binding Protein

(Cat. No. HIGFBMAG-53K)

IGFBP1	IGFBP5
IGFBP2	IGFBP6
IGFBP3°	IGFBP7
IGFBP4	

Human IGF

(Cat. No. HIGFMAG-52K)

IGF-1	IGF-2
-------	-------

Legend key for MILLIPLEX® kits

- † Analytes which cannot be plexed together:
 - Active and total
 - Free and total
- ◆◆ Available in Cat. No. listed
- ▼ Premix panel only
- Bulk indicates Space Saver Packaging
- ▲ These analytes cannot be plexed with other analytes in this panel in serum/plasma

- ‡ Requires a protease inhibitor during sample collection
- Requires sample extraction
- ** No standard/QC; antibody levels are reported as MFIs (Median Fluorescence Intensities)
- Serum/Plasma only
- ★ Tissue Culture samples only

- ^ Competitive assay format
- Available for custom premix
- ° Analyte requires a different sample dilution from others in panel
- Analyte cannot be measured in serum
- # Requires sample acidification
- ° Analyte cannot be measured in plasma
- u Urine only

Human Pituitary Panel 1

(Cat. No. HPTP1MAG-66K)

ACTH	FSH
Agouti-Related Protein (AgRP)	GH
CNTF	LH
	TSH

Cardiovascular**Human CVD Panel 1**

(Cat. No. HCVD1MAG-67K)

BNP	LIGHT
CK-MB	NT proBNP
CXCL6/GCP-2	Oncostatin-M (OSM)
CXCL16	Placental Growth Factor (PIGF)
Endocan-1 (ESM-1)	Troponin I (TnI)
FABP3	
FABP4	

Human CVD Panel 2

(Cat. No. HCVD2MAG-67K)

ADAMTS13	Myoglobin
D-dimer ^s	NGAL/Lipocalin-2
GDF15	sP-Selectin
sICAM-1	Serum Amyloid A
Myeloperoxidase (MPO)	sVCAM-1

Human CVD Panel 3 (Acute Phase)

(Cat. No. HCVD3MAG-67K)

α-1-Acid Glycoprotein (AGP)	Haptoglobin
α-2-Macroglobulin	sL-Selectin
Adipsin/Factor D	Platelet Factor 4 (PF4)/CXCL4
CRP	Serum Amyloid P (SAP)
Fetuin A	
Fibrinogen [⊙]	

Human CVD Panel 4

(Cat. No. HCVD4MAG-67K)

sCD31/sPECAM-1	Thrombomodulin
sE-Selectin	Tissue Factor (TF)
Follistatin (FST)	Troponin T (TnT)
Pentraxin-3 (PTX3)	

Human CVD Panel 6

(Cat. No. HCVD6MAG-67K)

sCD14	MCAM/MUC18
DPP4	Progranulin (pGRN)
Endostatin	ucMGP
LRG1	

Human Apolipoprotein

(Cat. No. APOMAG-62K)

Apo AI	Apo CII
Apo AII	Apo CIII
Apo B	Apo E

Bone**Human Bone**

(Cat. No. HBNMAG-51K)

ACTH	Osteocalcin (OC)
DKK1	Osteopontin (OPN)
FGF-23	Osteoprotegerin (OPG)
IL-1β	PTH
IL-6	Sclerostin (SOST)
Insulin	TNFα
Leptin	

Human RANKL – Singleplex

(Cat. No. HRNKLMAG-51K-01)

RANKL

Cancer**Human Immuno-Oncology Checkpoint Protein Panel 1**

(Cat. No. HCKP1-11K)

17 (Cat. No. HCKP1-11K-PX17)
 17 (Bulk Cat. No. HCKP1-11K-PXBK17)

BTLA	HVEM
CD27	ICOS
CD28	LAG-3
CD40	PD-1
CD80/B7-1	PD-L1
CD86/B7-2	PD-L2
CTLA-4/CD152	TIM-3
GITR	TLR-2
GITRL	

Human Immuno-Oncology Checkpoint Protein Panel 2

(Cat. No. HCKP2-11K)

4-1BBL/TNFSF9	FGL1/Hepassocin
5'-NT/CD73	Galectin-1
APRIL	Galectin-3
Arginase-1	Granulysin
B7-H2/ICOSL	Granzyme B
B7-H3/CD276	IDO1
B7-H4/VTCN1	MICA
B7-H5/VISTA	MICB
B7-H6	Nectin-2
BAFF/BLyS	Nectin-4
CD25/IL-2Rα	OX40/CD134
CD30/TNFRSF8	Perforin
CD40L	PVR/CD155
CD137/4-1BB	Siglec-7
CD226/DNAM-1	Siglec-9
E-Cadherin	

Human Cancer Autoantibody IgG**

(Cat. No. HCABMAG-13K)

QUALITATIVE KIT

CCNB1/Cyclin B1	IGF2BP2/IMP2
CTAG1B/NY-ESO-1	IGF2BP3/IMP3
ENO1	Mucin-1 (MUC1)
Galectin-1 (GAL1)	P16-INK4a
Galectin-3 (GAL3)	p53
HER2/ErbB2	SOX2
HIF-1α	Survivin/BIRC5
HSP60	

Human Exosome Characterization

(Cat. No. HEXSM-170K)

9 (Cat. No. HEXSM-170K-PMX)
 9 (Bulk Cat. No. HEXSM-170K-BKPMX)

Argonaute-2 (Ago2)*	Flotillin-1 (FLOT-1)
Calreticulin (CALR)	GAPDH
CD9	Syntenin-1 (Syn-1)
CD63	TSG101
CD81	

*The Ago2 bead is static and will always be included in the kit format.

Human Circulating Cancer Biomarker Panel 1

(Cat. No. HCCBP1MAG-58K)

α-Fetoprotein (AFP)	IL-8/CXCL8
CA125	Leptin
CA15-3	MIF
CA19-9	Osteopontin (OPN)
CEA	Prolactin
CYFRA21-1	PSA (free) [†]
sFAS/TNFRSF6	PSA (total) [†]
sFasL	SCF
FGF-2/FGF-basic	TGFα
HCGβ	TNFα
HE4	TRAIL/TNFSF10
HGF	VEGF-A
IL-6	

Human Circulating Cancer Biomarker Panel 3

(Cat. No. HCCBP3MAG-58K)

Cathepsin D	Melanoma Inhibitory Activity (MIA)
Ferritin	Myeloperoxidase (MPO)
Fibroblast Activation Protein (FAP)	Sex Hormone Binding Globulin (SHBG)
Galectin-3	
IGFBP3 [⊙]	

MILLIPLEX® Cancer Kits

See how our kits are advancing cancer and immuno-oncology research at
SigmaAldrich.com/milliplex-cancer



Human Circulating Cancer Biomarker Panel 4

(Cat. No. HCCB4MAG-58K)

ALDH1A1	Kallikrein-6
Carbonic Anhydrase 9 (CA9)	Mesothelin
CD44	Midkine
EpCAM	NCAM1/L1CAM/CD171
Hepsin	Transglutaminase 2 (TGM2)

Human Cancer/Metastasis Biomarker Panel 1

(Cat. No. HCMBMAG-22K)

DKK1	Osteoprotegerin (OPG)
GDF15	Periostin
Neuron-specific Enolase (NSE)	TRAP5
Osteonectin/SPARC	TWEAK
	YKL40/CHI3L1

Human Angiogenesis/Growth Factor Panel 1

(Cat. No. HAGP1MAG-12K)

Angiopoietin-2	HB-EGF
BMP-9	HGF
EGF	IL-8/CXCL8
Endoglin	Leptin
Endothelin-1	Placental Growth Factor (PlGF)
FGF-1/FGF-acidic	VEGF-A
FGF-2/FGF-basic	VEGF-C
Follistatin (FST)	VEGF-D
G-CSF	

Human Angiogenesis Panel 2

(Cat. No. HANG2MAG-12K)

Angiostatin/Kringle	sIL-6Rα
sAXL	sNeuropilin-1 (sNRP-1)
sCD31/sPECAM-1	Osteopontin (OPN)
sc-Kit/sStem Cell Factor Receptor (SCFR)	PDGF-AB/BB
sE-Selectin	Thrombospondin-2 (TSP-2)
sEGFR/sHER1/sErbB1	sTIE-2
sHER2/sEGFR2/sErbB2	suPAR
sHER3/sEGFR3/sErbB3	sVEGFR1/sFlt-1
sHGFR/sc-Met	sVEGFR2/sKDR/sFlk-1
	sVEGFR3/sFlt-4

Neuroscience**Human Amyloid Beta and Tau (CSF samples)**

(Cat. No. HNABTMAG-68K)

Amyloid beta 1-40	pTau (Thr181)
Amyloid beta 1-42	Tau (total)

Human Neuroscience Panel 1 (CSF samples)

(Cat. No. HNS1MAG-95K)

α-Synuclein	PARK5/UCHL1
Glial Fibrillary Acidic Protein (GFAP)	PARK7/DJ1
Neuron-specific Enolase (NSE)	Transglutaminase 2 (TGM2)

Human Neuroscience Panel 2

(Cat. No. HNS2MAG-95K)

Angiogenin (ANG) ^o	Ferritin
ApoE4 ^o	Neurogranin (NRGN)
FABP3	TREM2

Human Neurodegenerative Autoantibody IgG**

(Cat. No. HNDGAB-220K)

QUALITATIVE KIT

α-synuclein	DRD2
Aβ42	ELAVL4 (HuD)
AChR	GFAP
AGNA (SOX1)	MBP
Ago2 (Su)	mGluR1
AMPA-R	MOG
Amphiphysin	NF-L
ANNA2 (Ri)/NOVA1	NMDAR1
AQP4	NRXN3
C1q	PCA1 (Yo)/CDR2
CRMP5	PNMA1 (Ma1)
CRYAB	tTau
DPPX	

Human Neurodegenerative Disease Panel 1

(Cat. No. HNDG1MAG-36K)

α-2-Macroglobulin	Complement Factor C3
Apo AI	Complement Factor H (CFH)
Apo CIII	
Apo E	

Human Neurodegenerative Disease Panel 2

(Cat. No. HNDG2MAG-36K)

α-1-Antitrypsin (A1AT)	MIP-4/PARC/CCL18
C4	PEDF
CRP	Serum Amyloid P (SAP)

Human Neurodegenerative Disease Panel 3

(Cat. No. HNDG3MAG-36K)

BDNF	PAI-1 (total)
Cathepsin D	PDGF-AA
sICAM-1	PDGF-AB/BB
Myeloperoxidase (MPO)	RANTES/CCL5
sNCAM	sVCAM-1

Human Neurodegenerative Disease Panel 4 (CSF samples)

(Cat. No. HNDG4MAG-36K)

Amyloid beta 1-40	sRAGE
Amyloid beta 1-42	S100B
GDNF	

Human Neurological Disorders Panel 3

(Cat. No. HND3MAG-39K)

Angiotensinogen (AGT)	Osteopontin (OPN)
Contactin-1	Soluble Superoxide Dismutase 1 (sSOD1)
Fetuin A ^o	Soluble Superoxide Dismutase 2 (sSOD2)
Kallikrein-6	

Human Neuropeptide[■] ^

(Cat. No. HNPMAg-35K)

α-MSH	Oxytocin
β-Endorphin	Substance P
Neurotensin	

Human Circadian Stress[■] ^

(Cat. No. HNCMSAG-35K)

Cortisol	Melatonin
----------	-----------

Explore Biomarkers for Alzheimer's Disease Research

Learn how to combine multiplex and high sensitivity assay analysis for AD biomarker studies at SigmaAldrich.com/immunoassays-AD

**Legend key for MILLIPLEX® kits**

- † Analytes which cannot be plexed together:
 - Active and total
 - Free and total
- ♦ Available in Cat. No. listed
- ▼ Premix panel only
- Bulk indicates Space Saver Packaging
- ▲ These analytes cannot be plexed with other analytes in this panel in serum/plasma

- Requires a protease inhibitor during sample collection
- Requires sample extraction
- ** No standard/QC; antibody levels are reported as MFIs (Median Fluorescence Intensities)
- Serum/Plasma only
- ★ Tissue Culture samples only

- ^ Competitive assay format
- Available for custom premix
- ° Analyte requires a different sample dilution from others in panel
- Analyte cannot be measured in serum
- # Requires sample acidification
- ° Analyte cannot be measured in plasma
- u Urine only

Toxicity

Human Kidney Injury Panel 1 (Urine samples)

(Cat. No. HKI1MAG-99K)

Calbindin	KIM-1
Collagen IV	Renin
GSTa	TFF-3
IP-10/CXCL10	TIMP-1

Human Kidney Injury Panel 2 (Urine samples)

(Cat. No. HKI2MAG-99K)

α-1-Microglobulin	EGF
Albumin	NGAL/Lipocalin-2
Clusterin	Osteopontin (OPN)
Cystatin C	

Human Kidney Injury Panel 4 (Serum/Plasma samples)

(Cat. No. HKI4MAG-99K)

EGF	Osteopontin (OPN)
FABP1	PTH
IP-10/CXCL10	Renin
KIM-1	

Human Kidney Injury Panel 5 (Serum/Plasma samples)

(Cat. No. HKI5MAG-99K)

α-1-Microglobulin	Osteoactivin
Collagen IV	TIMP-1
NGAL/Lipocalin-2	Uromodulin

Human Kidney Injury Panel 6 (Serum/Plasma samples)

(Cat. No. HKI6MAG-99K)

β-2-Microglobulin	Cystatin C
Clusterin	RBP4

Human Liver Injury

(Cat. No. HLINJMAG-75K)

5'NT/CD73	MDH1
ARG1	SDH
GSTa	

PRIMATE

Immunology

Non-Human Primate Cytokine/ Chemokine /Growth Factor Panel A

☐ (Cat. No. PRCYTA-40K)

38 (Cat. No. PRCYTA-40K-PX38)♦

38 (Cat. No. PRCYTA-40K-BK38)♦

48 (Cat. No. PRCYTA-40K-PX48)

48 (Cat. No. PRCYTA-40K-BK48)

BCA-1/CXCL13	IL-15♦
sCD137/4-1BB/ TNFRSF9♦	IL-16
sCD40L♦	IL-17A/CTLA8♦
Eotaxin/CCL11	IL-17E/IL-25
sFASL♦	IL-18♦
FGF-2/FGF-basic	IL-21♦
Fractalkine/CX3CL1	IL-22♦
G-CSF♦	IL-23♦
GM-CSF♦	IL-28A/IFNλ2
Granzyme A♦	IL-31
Granzyme B♦	IL-33/NF-HEV (mature)♦
IFNα2♦	IP-10/CXCL10♦
IFNγ♦	I-TAC/CXCL11♦
IL-1α	MCP-1/CCL2♦
IL-1β♦	MIG/CXCL9♦
IL-1RA♦	MIP-1α/CCL3♦
IL-2♦	MIP-1β/CCL4♦
IL-4♦	MIP-3α/CCL20♦
IL-5♦	Perforin♦
IL-6♦	RANTES/CCL5♦▲°
IL-7♦	TGFα♦
IL-8/CXCL8♦	TNFα♦
IL-10♦	TNFβ/LTA
IL-12 (p70)♦	VEGF-A♦

Metabolism/ Endocrinology

Non-Human Primate Metabolic Hormone

(Cat. No. NHPMHMAG-45K)

Amylin (active)♦	Insulin
C-Peptide	Leptin
Ghrelin (active)♦	MCP-1/CCL2
GIP (total)	Pancreatic Polypeptide (PP)
GLP-1 (active)♦	PYY (total)
Glucagon♦	
IL-6	

Non-Human Primate Pituitary Panel 1

(Cat. No. NHPPT1MG-46K)

ACTH	FSH
Agouti-Related Protein (AgRP)	GH
CNTF	LH
	TSH

MOUSE

Immunology

Rat/Mouse Complement **NEW**

☐ (Cat. No. RMCMP-260K)

Adipsin/Factor D (Mouse)	Factor H (Rat, Mouse)
C1q (Mouse)	MBL1 (Rat, Mouse)
C1ra (Rat, Mouse)	MBL2 (Rat, Mouse)
C3b/iC3b+ (Mouse)	Properdin (Rat, Mouse)
C4b (Mouse)	Vitronectin (Mouse)
C5a (Rat, Mouse)	
Clusterin (Rat)	

Mouse Cytokine/Chemokine/ Growth Factor Expanded Panel 1

☐ (Cat. No. MCYT1-190K)

43 (Cat. No. MCYT1-190K-SPX)♦

43 (Cat. No. MCYT1-190K-SPXBK)♦

68 (Cat. No. MCYT1-190K-LPX)

68 (Cat. No. MCYT1-190K-LPXBK)

BAFF	IL-12p40♦
Betacellulin	IL-12p70♦
CCL25/TECK	IL-13♦
sCD137/4-1BB/ TNFRSF9	IL-15♦
CHI3L1	IL-16
CXCL16	IL-17A/CTLA8♦
EGF♦	IL-17F♦
Eotaxin/CCL11♦	IL-18♦
Erythropoietin (EPO)	IL-20
Exodus-2/ CCL21/6Ckine	IL-22♦
sFas/TNFRSF6	IL-28B/IFNλ3♦
sFasL	IL-31
FGF-2♦	IL-33♦
FLT3L	IP-10/CXCL10♦
Fractalkine/CX3CL1	KC/GRO/CINC-1/ CXCL1♦
G-CSF♦	LIF♦
GDF-15	LIX♦
GM-CSF♦	MDC/CCL22
Granzyme B	MCP-1/CCL2♦
sICAM-1	MCP-2
IFNα♦	MCP-3
IFNβ♦	MCP-5/CCL12
IFNγ♦	M-CSF♦
IL-1α♦	MIG/CXCL9♦
IL-1β♦	MIP-1α/CCL3♦
IL-2♦	MIP-1β/CCL4♦
IL-3♦	MIP-2/CXCL2♦
IL-4♦	MIP-3α/CCL20
IL-5♦	MIP-3β/CCL19
IL-6♦	sRAGE
IL-7♦	RANTES/CCL5♦
IL-9♦	TARC/CCL17
IL-10♦	TNFα♦
IL-11	VEGF-A♦

Non-Human Primate Cytokine Analysis

Save time and sample with the largest bead-based non-human primate multiplex assay on the market. Read more at [SigmaAldrich.com/milliplex-primate](https://www.sigmaaldrich.com/milliplex-primate).



Mouse High Sensitivity T Cell

- (Cat. No. MHSTCMAG-70K)
 18 (Cat. No. MHSTCMAG-70KPMX)
 18 (Bulk Cat. No. MHSTCMAG-70KPBK)

GM-CSF	IL-10
IFN γ	IL-12 (p70)
IL-1 α	IL-13
IL-1 β	IL-17A/CTLA8
IL-2	KC/GRO/CINC-1/ CXCL1
IL-4	LIX
IL-5	MCP-1/CCL2
IL-6	MIP-2/CXCL2
IL-7	TNF α

Mouse Th17

- (Cat. No. MTH17MAG-47K)
 25 (Cat. No. MT17MAG47K-PX25)
 25 (Bulk Cat. No. MT17MAG47PMX25BK)

sCD40L	IL-17E/IL-25
GM-CSF	IL-17F
IFN γ	IL-21
IL-1 β	IL-22
IL-2	IL-23
IL-4	IL-27
IL-5	IL-28B/IFN λ 3
IL-6	IL-31
IL-10	IL-33/NF-HEV (mature)
IL-12 (p70)	MIP-3 α /CCL20
IL-13	TNF α
IL-15	TNF β /Lymphotoxin- α (LTA)
IL-17A/CTLA8	

Mouse CD8+ T Cell

- (Cat. No. MCD8MAG-48K)
 15 (Cat. No. MCD8MAG48K-PX15)
 15 (Bulk Cat. No. MCD8MAG48KPBX15BK)

sCD137/4-1BB/ TNFRSF9	IL-4
sFas/TNFRSF6	IL-5
sFasL	IL-6
GM-CSF	IL-10
Granzyme B	IL-13
IFN γ	MIP-1 β /CCL4
IL-2	RANTES/CCL5
	TNF α

Mouse Soluble Cytokine Receptor

- (Cat. No. MSCRMAG-42K)

sCD30	sRAGE
sGP130	sTNF RI
sIL-1RI	sTNF RII
sIL-1RII	sVEGFR1/sFlt-1
sIL-2Ra	sVEGFR2/sKDR/sFlk-1
sIL-4R	sVEGFR3/sFlt-4
sIL-6R	

**Mouse MMP Panel 1
(Serum/Plasma samples)**

- (Cat. No. MMMP1MAG-79K)

MMP-2	MMP-8
MMP-3	

**Mouse MMP Panel 2
(Serum/Plasma samples)**

- (Cat. No. MMMP2MAG-79K)

proMMP-9	MMP-12
----------	--------

**Mouse MMP Panel 3
(Cell culture samples)**

- (Cat. No. MMMP3MAG-79K)

MMP-2	proMMP-9
MMP-3	MMP-12
MMP-8	

Immune Response**Mouse Immunoglobulin Isotyping**

- (Cat. No. MGAMMAG-300K)

IgA	IgG2b
IgG1	IgG3
IgG2a	IgM

**Mouse IgE – Singleplex
(Cat. No. MGAMMAG-300E)**

IgE

**Metabolism/
Endocrinology****Mouse Adipokine
(Serum/Plasma samples)**

- (Cat. No. MADKMAG-71K)

IL-6	PAI-1 (total)
Insulin	Resistin
Leptin	TNF α
MCP-1/CCL2	

**Mouse Adipocyte
(Cell culture samples)**

- (Cat. No. MADCYMAG-72K)

Adiponectin	PAI-1 (total)
IL-6	Resistin
Leptin	TNF α
MCP-1/CCL2	

**Mouse Adiponectin –
Singleplex**

- (Serum/Plasma samples)

- (Cat. No. MADPNMAG-70K-01)

Adiponectin

**Mouse Metabolic Hormone
Expanded**

- (Cat. No. MMHE-44K)

Amylin (active) $\blacktriangleright$	Leptin
C-Peptide 2	MCP-1/CCL2
Ghrelin (active) $\blacktriangleright$	Pancreatic Polypeptide (PP)
GIP (total)	PYY (total)
GLP-1 (active) $\blacktriangleright$	Resistin
GLP-1 (total) $\blacktriangleright$	Secretin
Glucagon $\blacktriangleright$	TNF α
IL-6	
Insulin	

Mouse Myokine

- (Cat. No. MMYOMAG-74K)

BDNF	Irisin
Erythropoietin (EPO) $\star$	LIF
FGF-21	Myostatin/GDF8
Follistatin-like Protein 1 (FSTL1)	Oncostatin-M (OSM)
IL-6	Osteocrin (OSTN)/ Musclin
IL-15	Osteonectin/SPARC

Mouse Aging Panel 1

- (Cat. No. MAGE1MAG-25K)

CTACK/CCL27	Jag1
FGF-21	Leptin
GnRH	Notch1
IL-6	TNF α

Mouse Pituitary

- (Cat. No. MPTMAG-49K)

ACTH	LH
BDNF	Prolactin
FSH	TSH
GH	

Mouse TSH – Singleplex

- (Cat. No: RMTSH-200K)

- (Bulk Cat. No. RMTSH-200K-10PK)

TSH

Legend key for MILLIPLEX® kits

- $\dagger$ Analytes which cannot be plexed together:
 – Active and total
 – Free and total
 $\blacklozenge$ Available in Cat. No. listed
 $\blacktriangledown$ Premix panel only
 Bulk indicates Space Saver Packaging
 $\blacktriangle$ These analytes cannot be plexed with other
 analytes in this panel in serum/plasma

- $\blacktriangleright$ Requires a protease inhibitor during sample
 collection
 $\blacksquare$ Requires sample extraction
 $\star\star$ No standard/QC; antibody levels are
 reported as MFIs (Median Fluorescence
 Intensities)
 $\bullet$ Serum/Plasma only
 $\star$ Tissue Culture samples only

- $\wedge$ Competitive assay format
 $\blacksquare$ Available for custom premix
 $\circ$ Analyte requires a different sample dilution
 from others in panel
 $\odot$ Analyte cannot be measured in serum
 $\#$ Requires sample acidification
 $\textcircled{S}$ Analyte cannot be measured in plasma
 $\textcircled{U}$ Urine only

Cardiovascular

Mouse CVD Panel 1

(Cat. No. MCVD1MAG-77K)

sCD31/sPECAM-1	ProMMP-9
sE-Selectin	sP-Selectin
sICAM-1	Thrombomodulin
PAI-1 (total)	

Mouse CVD Panel 2

(Cat. No. MCVD2MAG-77K)

sCD40L	Oncostatin M (OSM)
CXCL16	Placental Growth Factor 2
Endocan-1 (ESM-1)	Troponin I (TnI)
Follistatin (FST)	Troponin T (TnT)
LIGHT	

Mouse Acute Phase Panel 2

(Cat. No. MAP2MAG-76K)

α-1-Acid Glycoprotein (AGP)	CRP
α-2-Macroglobulin	Haptoglobin
Adipsin/Factor D	Serum Amyloid P (SAP)

Cancer

Humanized Mouse

☐ (Cat. No. HUMU-210K)

41 (Cat. No. HUMU-210K-PMX)

41 (Cat. No. HUMU-210K-PMXBK)

G-CSF (Human)	IL-12p40 (Human)
G-CSF (Mouse)	IL-12p40 (Mouse)
GM-CSF (Human)	IL-12p70 (Human)
GM-CSF (Mouse)	IL-12p70 (Mouse)
IFNγ (Human)	IL-15 (Human)
IFNγ (Mouse)	IL-17A/CTLA8 (Human)
IL-1α (Human)	IL-17A/CTLA8 (Mouse)
IL-1β (Human)	IL-18 (Human)
IL-1β (Mouse)	IL-18 (Mouse)
IL-2 (Human)	IP-10/CXCL10 (Human)
IL-2 (Mouse)	IP-10/CXCL10 (Mouse)
IL-3 (Human)	MCP-1/CCL2 (Human)
IL-3 (Mouse)	MCP-1/CCL2 (Mouse)
IL-4 (Human)	M-CSF (Human)
IL-5 (Human)	M-CSF (Mouse)
IL-5 (Mouse)	MIG/CXCL9 (Human)
IL-6 (Human)	MIG/CXCL9 (Mouse)
IL-6 (Mouse)	TNFA (Human)
IL-7 (Human)	TNFA (Mouse)
IL-8 (Human)	
KC/GRO/CINC-1/CXCL1 (Mouse)	
IL-10 (Human)	

Mouse Angiogenesis / Growth Factor Panel 1

(Cat. No. MAGPMAG-24K)

sALK-1	IL-6
Amphiregulin	IL-17A/CTLA8
Angiopoietin-2 ▲	KC/GRO/CINC-1/CXCL1
Betacellulin ▲	Leptin
sCD31/sPECAM-1 ▲	MCP-1/CCL2
EGF	MIP-1α/CCL3
Endoglin	Placental Growth Factor 2
Endothelin-1	Prolactin
sFasL	SDF-1/CXCL12
FGF-2/FGF-basic	TNFA
Follistatin (FST)	VEGF-A
G-CSF	VEGF-C
HGF	VEGF-D
IL-1β	

Mouse Immuno-Oncology Checkpoint Protein Panel 1

☐ (Cat. No. MCKP1-110K)

28 (Cat. No. MCKP1-110K-PX28)

28 (Bulk Cat. No. MCKP1-110K-PBK28)

4-1BBL/TNFSF9	Galectin-1
5'-NT/CD73	Galectin-3
B7-H2/ICOSL	GITR
B7-H3/CD276	Granzyme B
BCA-1/CXCL-13	HVEM
BTLA	IFNγ
CD137/4-1BB	IL-10
CD226/DNAM-1	LAG-3
CD25/IL-2Ra	PD-1
CD27	PD-L1
CD40	PD-L2
CD80/B7-1	TIM-3
CTLA-4/CD152	TLR-2
E-Cadherin	TNFA

Neuroscience

Mouse Amyloid Beta

(Cat. No. MABMAG-83K)

Amyloid beta 1-40	Amyloid beta 1-42
-------------------	-------------------

Mouse Neuropeptide ■ ^

(Cat. No. RMNPMAG-83K)

α-MSH	Oxytocin
β-Endorphin	Substance P
Neurotensin	

Toxicity

Mouse Kidney Injury Panel 1

(Cat. No. MKI1MAG-94K)

β-2-Microglobulin ^u	Renin
IP-10/CXCL10	TIMP-1
KIM-1	VEGF-A ^u

Mouse Kidney Injury Panel 2

(Cat. No. MKI2MAG-94K)

Clusterin	NGAL/Lipocalin-2
Cystatin C	Osteopontin (OPN)
EGF	

RAT

Immunology

Rat Cytokine/Chemokine

☐ (Cat. No. RECYTMAG-65K)

27 (Cat. No. RECYMAG65K27PMX)

27 (Bulk Cat. No. RECYMAG65PMX27BK)

EGF	IL-10
Eotaxin/CCL11	IL-12 (p70)
Fractalkine/CX3CL1	IL-13
G-CSF	IL-17A/CTLA8
GM-CSF	IL-18
KC/GRO/CINC-1/CXCL1	IP-10/CXCL10
IFNγ	Leptin
IL-1α	LIX
IL-1β	MCP-1/CCL2
IL-2	MIP-1α/CCL3
IL-4	MIP-2/CXCL2
IL-5	RANTES/CCL5
IL-6	TNFA
	VEGF

Metabolism/Endocrinology

Rat Adipokine (Serum/Plasma samples)

(Cat. No. RADPKMAG-80K)

IL-1β	MCP-1/CCL2
IL-6	PAI-1 (total)
Insulin	TNFA
Leptin	

Rat Adipocyte (Cell culture samples)

(Cat. No. RADPCMAG-82K)

Adiponectin	MCP-1/CCL2
IL-1β	PAI-1 (total)
IL-6	TNFA
Leptin	

Get Results You Can Trust

See why MILLIPLEX® assays set the standard in quality-driven multiplexing.

Read more at SigmaAldrich.com/milliplex-quality



Rat Metabolic Hormone Expanded

(Cat. No. RMHE-120K)

Amylin (active)†	Insulin
C-Peptide 2	Leptin
Ghrelin (active)†	MCP-1/CCL2
GIP (total)	Pancreatic Polypeptide (PP)
GLP-1 (active) +†	PYY (total)
GLP-1 (total) +†	Secretin
Glucagon†	TNFα
IL-6	

Rat Myokine

(Cat. No. RMYOMAG-88K)

BDNF	IL-15
Erythropoietin (EPO)	Irisin
FGF-21	LIF
Follistatin-like Protein 1 (FSTL1)	Myostatin/GDF8
Fractalkine/CX3CL1	Osteocrin (OSTN)/Musclin
IL-6	Osteonectin/SPARC

Rat Pituitary

(Cat. No. RPTMAG-86K)

ACTH	LH
BDNF	Prolactin
FSH	TSH
GH	

Rat Stress Hormone■

(Cat. No. RSHMAG-69K)

ACTH	Melatonin^
Corticosterone^	

Rat Thyroid

(Cat. No. RTHYMAG-30K)

T3^	TSH
T4^	

Rat TSH – Singleplex

(Cat. No. RMTSH-200K)

(Bulk Cat. No. RMTSH-200K-10PK)

TSH

Cardiovascular

Rat Cardiac Injury Panel 1

(Cat. No. RCI1MAG-87K)

Cardiac Troponin I (cTnI)	FABP3
Cardiac Troponin T (cTnT)	Follistatin-like Protein 1 (FSTL1)
Creatine Kinase Muscle (CKM)	Myosin Light Chain 3 (MYL3)
	TIMP-1

Rat Vascular Injury Panel 1 (Serum/Plasma samples)

(Cat. No. RV1MAG-26K)

Caveolin-1	IL-6
Connective Tissue Growth Factor (CTGF)	MCP-1/ CCL2
KC/GRO/CINC-1/CXCL1	PAI-1 (total)
	TIMP-1
	TNFα
	VEGF

Rat Vascular Injury Panel 2 (Serum/Plasma samples)

(Cat. No. RV2MAG-26K)

Adiponectin	sICAM-1
sE-Selectin	

Rat Vascular Injury Panel 3 (Serum/Plasma samples)

(Cat. No. RV3MAG-26K)

α-1-Acid Glycoprotein (AGP)	Haptoglobin
α-2-Macroglobulin (A2M)	

Bone

Rat Bone Panel 1 (Serum/Plasma samples)

(Cat. No. RBN1MAG-31K)

ACTH	Leptin
DKK1	Osteoprotegerin (OPG)
FGF-23	PTH
Insulin	Sclerostin (SOST)

Neuroscience

Rat Neuropeptide■^

(Cat. No. RMNPMAG-83K)

α-MSH	Oxytocin
β-Endorphin	Substance P
Neurotensin	

Toxicity

Rat Kidney Toxicity Panel 1 (Urine samples)

(Cat. No. RKT1MAG-37K)

Clusterin	Osteopontin (OPN)
GSTα^	TIMP-1
IP-10/CXCL10	VEGF-A
KIM-1	

Rat Kidney Toxicity Panel 2 (Urine samples)

(Cat. No. RKT2MAG-37K)

α-1-Acid Glycoprotein (AGP)	Cystatin C
Albumin^	EGF
β-2-Microglobulin	NGAL/ Lipocalin-2

Rat Liver Injury

(Cat. No. RLI1MAG-92K)

5'NT/CD73	GSTα
ARG1	SDH
GOT1	

Legend key for MILLIPLEX® kits

- + Analytes which cannot be plexed together:
 - Active and total
 - Free and total
- ◆ Available in Cat. No. listed
- ▼ Premix panel only
- Bulk indicates Space Saver Packaging
- ▲ These analytes cannot be plexed with other analytes in this panel in serum/plasma

- ▮ Requires a protease inhibitor during sample collection
- Requires sample extraction
- ** No standard/QC; antibody levels are reported as MFIs (Median Fluorescence Intensities)
- Serum/Plasma only
- ★ Tissue Culture samples only

- ^ Competitive assay format
- Available for custom premix
- Analyte requires a different sample dilution from others in panel
- Analyte cannot be measured in serum
- # Requires sample acidification
- ˆ Analyte cannot be measured in plasma
- u Urine only

COMPANION/ AGRICULTURAL

Immunology

Bovine Cytokine/Chemokine Panel 1

-  (Cat. No. BCYT1-33K)
 (Cat. No. BCYT1-33K-PX15)
 (Bulk Cat. No. BCYT1-33K-PXBK15)

IFN γ	IL-36RA/IL-1F5
IL-1 α	IP-10/CXCL10
IL-1 β	MCP-1/CCL2
IL-4	MIP-1 α /CCL3
IL-6	MIP-1 β /CCL4
IL-8/CXCL8	TNF α
IL-10	VEGF-A
IL-17A/CTLA8	

Bovine Acute Phase Protein

(Cat. No. BAPP1-150K)

Haptoglobin	Fibrinogen $\odot$
-------------	--------------------

Canine Cytokine/Chemokine

-  (Cat. No. CCYTOMAG-90K)
 (Cat. No. CCYTOMAG-90K-PX13)
 (Bulk Cat. No. CCYTOMAG90KPX13BK)

GM-CSF	IL-15
IFN γ	IL-18
IL-2	IP-10/CXCL10
IL-6	KC-like
IL-7	MCP-1/CCL2
IL-8/CXCL8	TNF α
IL-10	

Chicken Cytokine/Chemokine

-  (Cat. No. GCYT1-16K)
 (Cat. No. GCYT1-16K-PX12)
 (Cat. No. GCYT1-16K-PXBK12)

IFN α	IL-21
IFN γ	M-CSF
IL-2	MIP-1 β /CCL4
IL-6	MIP-3 α /CCL20
IL-10	RANTES/CCL5
IL-16	VEGF

Equine Cytokine/Chemokine

-  (Cat. No. EQCYTMAG-93K)
 (Cat. No. EQCYTMAG-93KPX23)
 (Cat. No. EQCTMG93KPX23BK)

Eotaxin/CCL11	IL-6
FGF-2/FGF-basic	IL-8/CXCL8
Fractalkine/CX3CL1	IL-10
G-CSF	IL-12 (p70)
GM-CSF	IL-13
GRO	IL-17A/CTLA8
IFN γ	IL-18
IL-1 α	IP-10/CXCL10
IL-1 β	MCP-1/CCL2
IL-2	RANTES/CCL5
IL-4	TNF α
IL-5	

Feline Cytokine/Chemokine ▼

-  (Cat. No. FCYTOMAG-20K-PMX)
 (Bulk Cat. No. FCYTOMAG20KPX19BK)
NON-CONFIGURABLE KIT

sFAS/TNFRSF6	IL-13
Flt3 Ligand	IL-18
GM-CSF	KC/GRO/CINC-1/ CXCL1
IFN γ	
IL-1 β	MCP-1/CCL2
IL-2	PDGF-BB
IL-4	RANTES/CCL5
IL-6	SCF
IL-8/CXCL8	SDF-1/CXCL12
IL-12 (p40)	TNF α

Ovine Cytokine/Chemokine Panel 1

-  (Cat. No. SCYT1-91K)
 (Cat. No. SCYT-91K-PX14)
 (Bulk Cat. No. SCYT-91K-PXBK14)

IFN γ	IL-17A/CTLA8
IL-1 α	IL-36RA/IL-1F5
IL-1 β	IP-10/CXCL10
IL-4	MIP-1 α /CCL3
IL-6	MIP-1 β /CCL4
IL-8/CXCL8	TNF α
IL-10	VEGF-A

Porcine Cytokine/Chemokine/ Growth Factor Expanded Panel **NEW**

-  (Cat. No. PIGCYT-250K)
 (Cat. No. PIGCYT-250K-SPX) ♦
 (Cat. No. PIGCYT-250K-SPBK) ♦
 (Cat. No. PIGCYT-250K-LPX)
 (Cat. No. PIGCYT-250K-LPBK)

CCL3L1	IL-10 ♦
CD31/PECAM-1	IL-12 ♦
CD34	IL-17A
EGF	IL-18 ♦
Eotaxin	IL-28B/IFN λ 3
GM-CSF ♦	IP-10/CXCL10
IFN α	MCP-1
IFN β	MIG/CXCL9
IFN γ ♦	MIP-1 β
IL-1 α ♦	MMP-1
IL-1 β ♦	MMP-9
IL-1RA ♦	PAI-1
IL-2 ♦	PDGF-BB
IL-4 ♦	TGF α
IL-6 ♦	TNF α ♦
IL-8/CXCL8 ♦	VEGF-A

Porcine Acute Phase Protein

(Cat. No. PAPP1-160K)

CRP	Haptoglobin
-----	-------------

Metabolism/ Endocrinology

Canine Gut Hormone

(Cat. No. CGTMAG-98K)

Amylin (total) ▶	Insulin
Ghrelin (active) ▶	Leptin
GIP (total)	Pancreatic
GLP-1 (active) ▶	Polypeptide (PP)
Glucagon	PYY (total)

Canine Pituitary Expanded Panel

(Cat. No. CANPIT-96K)

ACTH	LH
FSH	Prolactin
GH	TSH

MILLIPLEX® Veterinary Medicine and Animal Health Kits

Fetch the latest in animal health research with our growing selection of multiplex assays.

Discover more at SigmaAldrich.com/vetmed



Toxicity**Canine Kidney Toxicity
Expanded Panel 1
(Urine samples)**

(Cat. No. CKT1MAG-97K)

Clusterin	MCP-1/CCL2
Cystatin C	NGAL/Lipocalin-2
IL-8/CXCL8	Osteopontin (OPN)
KIM-1	

**Canine Kidney Toxicity
Panel 2
(Urine samples)**

(Cat. No. CKT2MAG-97K)

Albumin	RBP4
β-2-Microglobulin	TFF-3

MULTI-SPECIES**Immunology****Multi-Species TGFβ –
Singleplex #**

(Cat. No. TGFBMAG-64K-01)

(Bulk Cat. No. TGFBMAG-64K-01BK)

Species: C, E, F, H, M, NHP, Por, R, guinea pig, hamster, rabbit

TGFβ1

**Multi-Species TGFβ –
3 Plex▼#**

(Cat. No. TGFBMAG-64K-03)

Species: H, M, NHP, R, guinea pig, hamster, mini pig; additionally, TGF-β1 and TGF-β2 can be detected in the serum of C, E, and rabbit

NON-CONFIGURABLE KIT

TGFβ1	TGFβ3
TGFβ2	

**Metabolism/
Endocrinology****Multi-Species Hormone ■^**

(Cat. No. MSHMAG-21K)

Species: C, E, F, H, M, NHP, Por, R, guinea pig, hamster, rabbit

Cortisol	T3
Estradiol	T4
Progesterone	Testosterone

**KITS TO BE
PHASED OUT****Porcine Cytokine/Chemokine**

● (Cat. No. PCYTMAG-23K)

13 (Cat. No. PCYTMG-23K-13PX)

13 (Bulk Cat. No. PCYTMAG23PMX13BK)

GM-CSF	IL-6
IFNγ	IL-8/CXCL8
IL-1α	IL-10
IL-1β	IL-12
IL-1Ra	IL-18
IL-2	TNFa
IL-4	

**Get More Out of Your Assays!****Dual Reporter IgG/IgM Reagent Pack**

Cat. No. DRAABPK-240K

Use with Autoantibody Panels and the Vector-Borne Disease Panel.

Learn how to redefine your multiplex assays at
SigmaAldrich.com/milliplex-serology**Legend key for MILLIPLEX® kits**

† Analytes which cannot be plexed together:
– Active and total
– Free and total

♦ Available in Cat. No. listed

▼ Premix panel only

Bulk indicates Space Saver Packaging

▲ These analytes cannot be plexed with other analytes in this panel in serum/plasma

♦ Requires a protease inhibitor during sample collection

■ Requires sample extraction

** No standard/QC; antibody levels are reported as MFIs (Median Fluorescence Intensities)

● Serum/Plasma only

★ Tissue Culture samples only

^ Competitive assay format

● Available for custom premix

◦ Analyte requires a different sample dilution from others in panel

○ Analyte cannot be measured in serum

Requires sample acidification

s Analyte cannot be measured in plasma

u Urine only

Cell Signaling MILLIPLEX® Assays

Akt/mTOR (Phosphoprotein) – 11 Plex

(Cat. No. 48-611MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB	✓ (Ser473)	H, M, R
GSK3α	✓ (Ser21)	H, M, R
GSK3β	✓ (Ser9)	H, M, R
IGF1R	✓ (Tyr1135/1136)	H, M
IR	✓ (Tyr1162/1163)	H
IRS1	✓ (Ser636)	H, M, R
mTOR	✓ (Ser2448)	H, M, R
p70S6 Kinase	✓ (Thr412)	H, M, R
PTEN	✓ (Ser380)	H, M, R
RPS6	✓ (Ser235/236)	H, M, R
TSC2	✓ (Ser939)	H, M, R

Akt/mTOR (Total) – 11 Plex

(Cat. No. 48-612MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB	✓	H, M, R
GSK3α	✓	H, M, R
GSK3β	✓	H, M, R
IGF1R	✓	H, M, R
IR	✓	H, R
IRS1	✓	H, M, R
mTOR	✓	H, M, R
p70S6 Kinase	✓	H, M, R
PTEN	✓	H, M, R
RPS6	✓	H, M, R
TSC2	✓	H, M, R

Early Apoptosis – 7 Plex

(Cat. No. 48-669MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB	✓ (Ser473)	H, M, R
BAD	✓ (Ser112)	H, M
Bcl-2	✓ (Ser70)	H
Active Caspase-8 (Asp384)	✓	H
Active Caspase-9 (Asp315)	✓	H
JNK/SAPK1	✓ (Thr183/Tyr185)	H, M, R
p53	✓ (Ser46)	H

Human Bcl-2 Family Apoptosis Panel 1 - 6 Plex

(Cat. No. 48-682MAG) AB1 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
BAD	✓	✓ (Ser112) H
BAX	✓	H
Bcl-xL/BAD (interaction)	✓	H
BIM	✓	H
Mcl-1/BIM (interaction)	✓	H

Human Bcl-2 Family Apoptosis Panel 2 - 4 Plex

(Cat. No. 48-683MAG) AB1 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Bcl-2	✓	H
Bcl-xL	✓	H
Mcl-1	✓	H
NOXA/Mcl-1 (interaction)	✓	H

Human Cell Health - 16 Plex

(Cat. No. 48-691MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
CHOP	✓	H
Cleaved PARP	✓	H
Complex IV	✓	H,M,R
Cyclin B1	✓	H
eIF4B		✓ (Ser422) H,M,R
GRP78	✓	H,M,R
HIF-1α	✓	H
Histone H2A.x		✓ (Ser139) H
Histone H3		✓ (Ser10) H
HSP70	✓	H
Ki-67	✓	H
LC3B	✓	H,M,R
NF-κB		✓ (Ser536) H
p53		✓ (Ser15) H
Polyubiquitin K48-linkage	✓	H,M
RNA Pol II		✓ (Ser2) H,M,R

Learn More About Cancer and Metastasis Signaling Pathways

Read about all of the kits that can elucidate your signaling research at SigmaAldrich.com/milliplex-cancer-cs



- Can be plexed with other 2 Plexes
- * Cannot plex with other phospho Akt
- ▼ Premix panel only

AB1: Uses Assay Buffer 1
AB2: Uses Assay Buffer 2

H Human
M Mouse
R Rat

Human DNA Damage/Genotoxicity – 7 Plex(Cat. No. 48-621MAG) AB1 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
ATR	✓	H
Chk1		✓ (Ser345) H, M, R
Chk2		✓ (Thr68) H
H2A.X		✓ (Ser139) H, M, R
MDM2	✓	H, M, R
p21	✓	H
p53		✓ (Ser15) H

MAPK/SAPK (Phosphoprotein) – 10 Plex(Cat. No. 48-660MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
ATF2		✓ (Thr71) H, M
c-Jun		✓ (Ser73) H, M, R
Erk/MAPK 1/2		✓ (Thr185/Tyr187) H, M, R
HSP27		✓ (Ser78) H
JNK/SAPK1		✓ (Thr183/Tyr185) H, M, R
MEK1		✓ (Ser222) H, M, R
MSK1		✓ (Ser212) H, M
p38/SAPK2A/B		✓ (Thr180/Tyr182) H, M, R
p53		✓ (Ser15) H
STAT1		✓ (Tyr701) H, M

Human RTK (Phosphoprotein)

(Choose analytes that meet your needs)

(Cat. No. HPRTKMAG-01K) AB1

Analyte	Total	Phosphorylated
c-Kit		(pan Tyr) H
c-Met/HGFR		(pan Tyr) H
ErbB1/EGFR		(pan Tyr) H
ErbB2/HER2		(pan Tyr) H
ErbB3/HER3		(pan Tyr) H
ErbB4/HER4		(pan Tyr) H
FGFR1		(pan Tyr) H
Flt3		(pan Tyr) H
IGF1R		(pan Tyr) H
M-CSFR		(pan Tyr) H
PDGFRβ		(pan Tyr) H
TIE1		(pan Tyr) H
TIE2		(pan Tyr) H
VEGFR1/Flt-1		(pan Tyr) H
VEGFR2/KDR/Flk-1		(pan Tyr) H
VEGFR3/Flt-4		(pan Tyr) H, M, R

Multi-Pathway (Phosphoprotein) – 9 Plex(Cat. No. 48-680MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB		✓ (Ser473) H, M, R
CREB		✓ (Ser133) H, M, R
Erk/MAPK 1/2		✓ (Thr185/Tyr187) H, M, R
JNK/SAPK1		✓ (Thr183/Tyr185) H, M, R
NFκB		✓ (Ser536) H, M, R
p38/SAPK2A/B		✓ (Thr180/Tyr182) H, M, R
p70S6 Kinase		✓ (Thr412) H, M, R
STAT3		✓ (Ser727) H, M, R
STAT5A/B		✓ (Tyr694/699) H, M, R

Multi-Pathway (Total) – 9 Plex(Cat. No. 48-681MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB	✓	H, M, R
CREB	✓	H, M, R
Erk/MAPK 1/2	✓	H, M, R
JNK/SAPK1	✓	H, M, R
NFκB	✓	H, M, R
p38/SAPK2A/B	✓	H, M, R
p70S6 Kinase	✓	H, M, R
STAT3	✓	H, M, R
STAT5A/B	✓	H, M, R

NFκB – 6 Plex(Cat. No. 48-630MAG) AB1 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
c-Myc	✓	H
FADD		✓ (Ser194) H
IκBα		✓ (Ser32) H
IKKα/β		✓ (Ser177/181) H
NFκB		✓ (Ser536) H, M
TNFR1	✓	H

Protein Translation – 6 Plex(Cat. No. 48-655MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
eIF2α		✓ (Ser51) H, M, R
eIF-4B		✓ (Ser422) H, M, R
eIF-4E		✓ (Ser209) H, M, R
eIF-4G		✓ (Ser1108) H, M, R
4E-BP1	✓	✓ (Thr37/46) H, M, R

Master Your Multiplexing

Get MILLIPLEX® multiplexing tips & tricks from the experts at SigmaAldrich.com/milliplex-tips



Human Ras-Raf Oncoprotein — 6 Plex(Cat. No. 48-684MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
BRAF	✓ (Ser446)	H, M, R
CRAF	✓ (Ser338)	H, M, R
MEK1	✓ (Ser217/221)	H, M, R
Ras	✓	H, M, R
RasG12D	✓	H
RasG12V	✓	H

Src Family Kinase Active Site (Phosphoprotein) — 8 Plex(Cat. No. 48-650MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Blk	✓ (Tyr389)	H
Fgr	✓ (Tyr412)	H, R
Fyn	✓ (Tyr420)	H, M, R
Hck	✓ (Tyr411)	H, M
Lck	✓ (Tyr394)	H, M, R
Lyn	✓ (Tyr397)	H, R
Src	✓ (Tyr419)	H, M, R
Yes	✓ (Tyr421)	H, M, R

STAT (Phosphoprotein) — 5 Plex(Cat. No. 48-610MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
STAT1	✓ (Tyr701)	H, M
STAT2	✓ (Tyr690)	H
STAT3	✓ (Tyr705)	H, M, R
STAT5A/B	✓ (Tyr694/699)	H, M, R
STAT6	✓ (Tyr641)	H

T-Cell Receptor (Phosphoprotein) — 7 Plex(Cat. No. 48-690MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
CD3ε	✓ (pan Tyr)	H
CREB	✓ (Ser133)	H, M, R
Erk/MAPK 1/2	✓ (Thr185/Tyr187)	H, M, R
LAT	✓ (pan Tyr)	H
Lck	✓ (pan Tyr)	H, M, R
Syk	✓ (pan Tyr)	H
ZAP-70	✓ (pan Tyr)	H

TGFβ — 6 Plex(Cat. No. 48-614MAG) AB2 ▼ **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB	✓ (Ser473)	H, M, R
Erk/MAPK 1/2	✓ (Thr185/Tyr187)	H, M, R
SMAD2	✓ (Ser465/467)	H, M, R
SMAD3	✓ (Ser423/425)	H, M, R
SMAD4	✓	H, M, R
TGFβRII	✓	H

Akt Phospho/Total — 2 Plex(Cat. No. 48-618MAG) AB2 ▼ * ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt/PKB	✓ (Ser473)	H, M, R
Akt/PKB	✓	H, M, R

Akt1 Phospho/Total — 2 Plex(Cat. No. 48-631MAG) AB2 ▼ * ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt1	✓ (Ser473)	H, M, R
Akt1	✓	H, M, R

Akt2 Phospho/Total — 2 Plex(Cat. No. 48-632MAG) AB2 ▼ * ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Akt2	✓ (Ser474)	H, M, R
Akt2	✓	H, M, R

CREB Phospho/Total — 2 Plex(Cat. No. 48-628MAG) AB2 ▼ ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
CREB	✓ (Ser133)	H, M, R
CREB	✓	H, M, R

Erk/MAPK 1/2 Phospho/Total — 2 Plex(Cat. No. 48-619MAG) AB2 ▼ ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated
Erk/MAPK 1/2	✓ (Thr185/Tyr187)	H, M, R
Erk/MAPK 1/2	✓	H, M, R

- Can be plexed with other 2 Plexes
- * Cannot plex with other phospho Akt
- ▼ Premix panel only

AB1: Uses Assay Buffer 1
AB2: Uses Assay Buffer 2

H Human
M Mouse
R Rat

JNK Phospho/Total – 2 Plex(Cat. No. 48-622MAG) AB2 ▼ ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated	
JNK/SAPK1		✓ (Thr183/Tyr185)	H, M, R
JNK/SAPK1	✓		H, M, R

mTOR Phospho/Total – 2 Plex(Cat. No. 48-625MAG) AB2 ▼ ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated	
mTOR		✓ (Ser2448)	H, M, R
mTOR	✓		H, M

p38 Phospho/Total – 2 Plex(Cat. No. 48-624MAG) AB1 or AB2 ▼ ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated	
p38/SAPK2A/B		✓ (Thr180/Tyr182)	H, M, R
p38/SAPK2A/B	✓		H, M, R

STAT3 Phospho/Total – 2 Plex(Cat. No. 48-623MAG) AB2 ▼ ●● **NON-CONFIGURABLE KIT**

Analyte	Total	Phosphorylated	
STAT3		✓ (Tyr705)	H, M, R
STAT3	✓		H, M, R

MAPmate™ Singleplex Loading Control Assays

Plex loading control MAPmate™ assays into existing MILLIPLEX® Cell Signaling Panels. Our Cell Signaling Buffer and Detection Kit (Cat. No. 48-602MAG), may be purchased if additional buffers are required.

MAPmate™ Kits	Species Homology (Buffer)	Cat. No.
β-Tubulin (total)	H, M, R (AB1 & AB2)	46-713MAG
GAPDH (total)	H (AB2)	46-667MAG

What's In A Typical Qualified Cell Signaling MILLIPLEX® Kit?

- Can be plexed with other 2 Plexes
- * Cannot plex with other phospho Akt
- ▼ Premix panel only

AB1: Uses Assay Buffer 1
AB2: Uses Assay Buffer 2

H Human
M Mouse
R Rat

Luminex® Multiplexing Platforms

Driving Assay Consistency

Combined with the largest portfolio of multiplex analytes available, we provide you the power of Luminex® xMAP® technology, the most trusted, widely used platform for biomarker screening and protein analysis. As a Luminex® partner, we are a preferred distributor of Luminex® instruments, accessories, and software.

Luminex® xMAP® technology offers many advantages compared to other immunoassay methods:

- **Multiplexing:** Data generated simultaneously is less prone to user error, enabling the end-user to construct cleaner data sets to explore multiple parameters.
- **Low sample volume:** With minimal hands-on time, screen multiple analytes in a single sample using a maximum of 25 µL per well.
- **Magnetic bead-based format:** Responds rapidly and efficiently to a magnetic field, enabling better and faster washing techniques, including high-throughput washing options.
- **Reproducibility:** High-volume production of xMAP® microspheres allows assay standardization that solid-phased flat arrays cannot provide.
- **Speed/high-throughput:** Simultaneously measure the concentration of a large number of different analytes in a single sample, in 96- or 384-well format, to do your work faster, gaining early and comprehensive data.



Specifications

Instruments	xMAP® INTELLIFLEX® DRSE*	xMAP® INTELLIFLEX®	FLEXMAP 3D®	Luminex® 200™
Cat. No.	INTELLIFLEX-DRSE	INTELLIFLEX-RU0	40-014	40-012
Multiplexing	Up to 500	Up to 500	Up to 500	Up to 100 (80 on MagPlex® microspheres)
Dynamic Range**	≥5.5 logs (RP1) ≥4.5 logs (RP2)	≥5.5 logs (RP1)	≥4.5 logs (RP1)	≥3.5 logs (RP1)
Read Time	96-well in ~20 min 384-well in ~75 min	96-well in ~20 min 384-well in ~75 min	96-well in ~20 min 384-well in ~75 min	96-well in ~45 min
Automation-Compatible Hardware	Front and side eject	Front eject only	Front eject only	Front eject only
Reporter Laser	532 nm (green) 405 nm (violet)	532 nm (green)	532 nm (green)	532 nm (green)
Plate Format	96-well and 384-well	96-well and 384-well	96-well and 384-well	96-well
Footprint (Linear Bench Space, including PC)	58.4 cm x 61 cm x 76.2 cm	58.4 cm x 61 cm x 76.2 cm	110 cm x 62 cm x 63 cm	115 cm x 60 cm x 50 cm
Weight	54.4 kg (120 lbs)	54.4 kg (120 lbs)	77.1 kg (170 lbs)	49 kg (113 lbs)

*DRSE = dual reporter side eject; **RP = reporter

Luminex® 200™ and FLEXMAP 3D® systems run on the Luminex® xPONENT® 4.3 data acquisition platform; the xMAP® INTELLIFLEX® systems run on INTELLIFLEX 1.0. All systems provide a consistent experience, from data input to results, no matter which instrument is used. All software adapts to fit the capabilities of each platform.

Next Generation Multiplexing


Harness the power of the xMAP® INTELLIFLEX® system at SigmaAldrich.com/intelliflex



Conferma® ELISAs

Conferma® ELISAs are developed with a focus on critical reagents to deliver assays with excellent detectability, sensitivity, and lot-to-lot consistency.

Analyte	Std. Range (pg/mL)	LLOQ (Curve) (pg/mL)	Cat. No.
IL-6 (Human)	1.17 – 150	1.17	EZIL6-98K
IL-8 (Human)	0.82 – 200	0.82	EZHIL8-100K
MCP-1 (Human)	7.81 – 1,000	7.81	EZMCP1-99KRM
TNF-α (Human)	1.65 – 400	1.65	EZHTNFA-150K



Three lots of the TNF-α Conferma® ELISAs were evaluated for parallelism. Variation was minimal demonstrating consistent lot-to-lot performance.

Curve	Parallelism
Lot 2 vs. Lot 1	0.983
Lot 3 vs. Lot 1	0.992

ELISAs Designed For Consistency

Get details on Conferma® ELISA performance and see the difference for yourself at: SigmaAldrich.com/conferma



Millipore® Brand ELISAs

Metabolic/Endocrine/ ELISAs	Species	Standard Curve Range	Sensitivity	Sample Volume	Cat. No.	Bulk Packaging Cat. No.*
Adiponectin	Human	1.56–200 ng/mL	0.2 ng/mL	10 µL	EZHADP-61K	EZHADP-61BK
Adiponectin (High Molecular Weight)	Human	1.5–200 ng/mL	0.5 ng/mL	10 µL	EZHMWAN-65K	
Adiponectin	Mouse	1–50 ng/mL	0.2 ng/mL	10 µL	EZMADP-60K	EZMADP-60BK
Amylin (active)	Human	1–100 pM	0.7 pM	50 µL	EZHA-52K	
C-Peptide	Human	0.2–20 ng/mL	0.05 ng/mL	10 µL	EZHCP-20K	EZHCP-20BK
C-Peptide	Canine	0.2–10 ng/mL	0.24 ng/mL	25 µL	EZCCP-47K	
C-Peptide 2	Mouse, Rat	25–1,600 pM	15.0 pM	20 µL	EZRMCP2-21K	
FGF-21	Human	31.25–2,000 pg/mL	10.0 pg/mL	50 µL	EZHFGF21-19K	
FGF-21	Mouse, Rat	49.4–12,000 pg/mL	10.0 pg/mL	10 µL	EZRMFGF21-26K	
FGF-23	Human	9.9–2,400 pg/mL	3.5 pg/mL	50 µL	EZHFGF23-32K	EZHFGF23-32BK
Ghrelin (active)	Human	25–2,000 pg/mL	15.0 pg/mL	20 µL	EZGRA-88K	EZGRA-88BK
Ghrelin (active)	Mouse, Rat	25–2,000 pg/mL	8.0 pg/mL	20 µL	EZRGRA-90K	
Ghrelin (total)	Human	100–5,000 pg/mL	50.0 pg/mL	20 µL	EZGRT-89K	EZGRT-89BK
Ghrelin (total)	Mouse, Rat	0.1–10 ng/mL	0.04 ng/mL	20 µL	EZRGRT-91K	
GIP (total)	Human	8.2–2,000 pg/mL	4.2 pg/mL	20 µL	EZHGIP-54K	EZHGIP-54BK
GIP (total)	Mouse, Rat	8.2–2,000 pg/mL	4.2 pg/mL	10 µL	EZRMGIP-55K	EZRMGIP-55BK
GLP-1 (active)	Multi-Species	2–100 pM	1.0 pM	100 µL	EGLP-35K	EGLP-35BK
GLP-1 High Sensitivity (active) Δ●●	Multi-Species	See data sheet	0.14 pM	50 µL	EZGLPHS-35K	EZGLPHS-35BK
GLP-1 (total)	Multi-Species	4.1–1,000 pM	1.0 pM	20–50 µL	EZGLP1T-36K	EZGLP1T-36BK
GLP-2	Multi-Species	1–64 ng/mL	0.3 ng/mL	50 µL	EZGLP2-37K	

Metabolic/Endocrine/ ELISAs	Species	Standard Curve Range	Sensitivity	Sample Volume	Cat. No.	Bulk Packaging Cat. No.*
Glucagon Δ	Human, Mouse, Rat	0.02–2 ng/mL	0.003 ng/mL	150–300 μ L	EZGLU-30K	EZGLU-30BK
Growth Hormone (GH)	Mouse, Rat	0.7–50 ng/mL	0.07 ng/mL	10 μ L	EZRMGH-45K	EZRMGH-45BK
Insulin	Human	2–200 μ U/mL	1.0 μ U/mL	20 μ L	EZHI-14K	EZHI-14BK
Insulin	Mouse, Rat	0.2–10 ng/mL	0.1 ng/mL	10 μ L	EZRMI-13K	EZRMI-13BK
Leptin	Canine	0.78–50 ng/mL	0.21 ng/mL	20 μ L	EZCL-31K	
Leptin	Mouse	0.2–30 ng/mL	0.05 ng/mL	10 μ L	EZML-82K	EZML-82BK
Leptin	Rat	0.2–30 ng/mL	0.08 ng/mL	10 μ L	EZRL-83K	EZRL-83BK
Leptin “Dual Range”	Human	0.5–100 ng/mL	0.2 ng/mL	25 μ L	EZHL-80SK	EZHL-80BK
Procollagen Type IIA N-Propeptide (PIIANP)	Human	lot dependent	30.0 ng/mL	5 μ L	EZPIIANP-53K	
Proinsulin (total)	Human	2–200 pM	0.5 pM	20 μ L	EZHPI-15K	EZHPI-15BK
PYY (total)	Human	10–2,000 pg/mL	6.5 pg/mL	20 μ L	EZHPPYT66K	
Resistin	Human	0.16–5 ng/mL	0.02 ng/mL	20 μ L	EZHR-95K	EZHR-95BK
SAA-3	Mouse	0.078–5 μ g/mL	0.078 μ g/mL	10 μ L	EZMSAA3-12K	
Neuroscience ELISAs	Species	Standard Curve Range	Sensitivity	Sample Volume	Cat. No.	
BDNF (Brain-Derived Neurotrophic Factor)	Human, Rat	7.8–500 pg/mL	7.8 pg/mL	50 μ L	CYT306	
NPY (Neuropeptide Y)	Human	5–1,000 pg/mL	2 pg/mL	50 μ L	EZHNPY-25K	
NPY (Neuropeptide Y)	Mouse, Rat	0.01–2 ng/mL	0.004 ng/ mL	20 μ L	EZRMNPY-27K	
PEDF (Pigment Epithelium-Derived Factor)	Human	0.9–62.5 ng/mL	0.9 ng/mL	50 μ L	CYT420	
Phosphorylated Neurofilament, (pNF-H) Sandwich	Multi- Species	0.0293–15 ng/mL	0.0585 ng/ mL	1–10 μ L	NS170	
S100B	Human	2.7–2,000 pg/mL	1.3 pg/mL	50 μ L	EZHS100B-33K	

* Bulk packaging available on select kits — more environmentally friendly and saves space (10 kit equivalent)

•• Preferred assay for measuring GLP-1 (active)

Δ Chemiluminescent assay

Sigma-Aldrich® Brand ELISAs

We offer a portfolio of >1,300 ELISAs, covering a diverse range of species, research areas, and sample types. Find your ELISA at SigmaAldrich.com/elisa



Ancillary Equipment to Optimize Your Assays



Instruments	BioTek® 405™ TSUV Plate Washer for SMC®/MILLIPLEX® Technology	BioTek® 405™ TS Magnetic 96-Well Washer
Cat. No.	95-0004-06	40-092AB
Supported Technology	SMC® immunoassays MILLIPLEX® immunoassays	MILLIPLEX® immunoassays
Features	• Custom magnetic plate specific for use with SMC® technology • Ultrasonic Advantage™ program and 4 buffer switching • Optimized cell wash including Shaking & Soaking	• Supports 96-well magnetic bead washing • Ultrasonic Advantage™ program for maintenance • Quick menu for priming and washing
Wash Protocols	Pre-programmed with SMC® and MILLIPLEX® assay wash protocols	Predefined sample methods for easy operation
Plate Format	96-well for SMC® assays 96- and 384-well for MILLIPLEX® assays	96-well plates
Washer/Dispenser Options	Vacuum filtration, biomagnetic separation	Biomagnetic separation
Volume Range	50 – 3,000 µL/well	25 – 3,000 µL/well, in 1 µL increments
Residual Volume	≤ 2 µL per well	≤ 2 µL per well
Wash or Read Speed	300 µL/well, 3 cycles; ≤ 30 s	300 µL/well, 3 cycles; ≤ 30 s
Wash Cycles	1 to 250	1 to 250
Shaking	Programmable up to 60 min: slow, medium, fast, or variable	Programmable up to 60 min: slow, medium, fast, or variable
Automation Capability	BioSpa, BioStack, BenchCel, 3rd Party	BioTek® BioStack microplate stacker and 3rd Party
Footprint	14" W x 17" D x 10" H (35.6 x 43.2 x 25.4 cm)	14" W x 17" D x 10" H (35.6 x 43.2 x 25.4 cm)
Weight	With internal buffer switching: 36 lb (16.5 kg)	30 lb (13.5 kg)

Related Products

Description	Cat. No.
BioTek® 50™ TS Microplate Washer	40-301
BioTek® 800™ TS Absorbance Reader with Gen5™ Software	40-006
BioTek® MultiFlo™ FX Multimode Dispenser	40-999
SMC® Plate Washer Evaluation Kit	03-0165-00
Flat Magnet, for 96-Well Plates, for use with MILLIPLEX® Assays (and for use with Cat. No. 95-0004-06)	7103016
Flat Magnet, for 384-Well Plates, for use with MILLIPLEX® Assays (and for use with Cat. No. 95-0004-06)	7103017
Handheld Magnetic Bead Washer Separator Block for 96-Well Flat Bottom or Conical Well Plates	40-285
BioTek® 800™ TS Product Qualification Package (to be sold with Cat. No. 40-006)	1560516
Absorbance Test Plate (to be sold with Cat. No. 40-006)	7260522
Vacuum Filter 96-Well (to be sold with Cat. No. 40-092AB)	1170008

Enhance Your Workflows

Explore BioTek® washers for our platforms at
SigmaAldrich.com/smc-washer and SigmaAldrich.com/milliplex-washer



Belysa® Immunoassay Curve Fitting Software

(Cat. No. 40-122)

Our proprietary software package for curve fitting, examination, and comparison of data acquired from a variety of immunoassay platforms. Belysa® software has been designed to provide an easy-to-use tool to optimize curve fit and confirm the integrity of data within a plate or across plates from multiple runs.

Key Features and Benefits


- User-friendly interface
- Automatic and manual curve fitting (4pl, 5pl, linear, competitive, cubic spline)
- Curve comparison tools
- Rules-based data flagging
- Suitable for use with RUO Luminex® instruments, FemtoQuest® platform, and ELISA .csv acquisition files
- Windows 10 and macOS (beta) compatible

Data displayed after selecting curve fit model:


- Sample concentration
- Standard recovery
- Equation of line
- Limit of Detection (LOD)
- Lower Limit of Quantitation (LLOQ)
- Minimum Detectable Dose (MD)

Other features:

- Data is exportable to: Excel, .csv, .txt, and .pdf
- Data autoflagging based on bead count, %CV, or BLOQ
- Curve comparison tool to evaluate parallelism and relative potency of standards



Example 1. Simply drag and drop the .csv output file into the Belysa® browser for instant data review.



Example 2. Comparison of two standard curves in Belysa® software from two lots of the same MILLIPLEX® assay run 1 year apart.

Ready To Start Comparing Curves?

Get more details about Belysa® Software at: SigmaAldrich.com/belysa



AAWmazing Automation

Automate the assays you already run with plug-and-play automation



Whether you run **MILLIPLEX®** multiplex panels, **Conferma®** ELISAs, or **SMC®** assays for the **FemtoQuest®** ultrasensitive system, the **AAW™ Automated Assay Workstation** automates your existing workflow with no protocol overhaul required.

Beginner-ready automation powered by **Opentrons®**. Flexibility, consistency, and efficiency all in one workstation.

Works With What You Have

Compatible with **MILLIPLEX®**, **ELISA**, and **SMC®** assays. Gain access to the **Millipore® Protocol Library** for verified ready-to-go automated assay protocols.

Rely On Your Results

Reduce manual variability across operators and sites. Supports multi-site comparability and reproducibility in long-term studies.

Maximize Walkaway Time

Run 96-well workflows with reduced hands-on time and consistent, reproducible results — plate-to-plate, day-to-day.

Modular and Affordable

Easily swap pipettes yourself or add modules as your workflow grows. Simple setup, extensive tech support, and cost-effective pricing.

Highlighted Compatible Assays



MILLIPLEX® Multiplex Assays



Conferma® ELISAs



SMC® Ultrasensitive Assays



Duolink® PLA

Our Assay Ready Workstation (Cat. No. 99900256) comes with the following:


- AAW™ Automated Assay Workstation
- AAW™ 1-Channel Pipette, 1,000 µL
- AAW™ 8-Channel Pipette, 1,000 µL
- AAW™ Flex Gripper
- Opentrons Flex™ Heater-Shaker (includes all adapters)
- Opentrons Flex™ Magnetic Block
- Opentrons Flex™ Deck Expansion
- Access to the **Millipore® Protocol Library** of verified assay protocols

These are all the accessories needed to run the verified protocols for **SMC®** ultrasensitive immunoassays, **MILLIPLEX®** multiplex immunoassays, and **Duolink® PLA** assays that are in the **Millipore® Protocol Library**. Kit protocols will continue to grow. Check our website for current list.

Get Automated Today

Find more protocols and information at SigmaAldrich.com/AAW-Yeah





Beyond Biomarkers

Keep your focus on the next stage

Rigorous Assay Development

Specificity ●●

- Antibody specificity is tested to ensure negligible cross-reactivity in the tested sample types
- Assay specificity is ensured by confirming consistent performance in singleplex and multiplex formats

Method Comparisons ●●

- We test our kits against in-house kits and kits from other vendors when available, ensuring the highest performing kits

Selectivity ●●

- Assay buffers and bead diluents are optimized to enhance antibody binding to only those analytes of interest in samples

Detection and Sensitivity ●

- The assay range and minimum detectable concentrations (minDCs or LLOQ) for each analyte is published in our protocols, and is based on actual samples

Performance in a Sample Matrix ●

- Assays are verified for the sample type noted in the protocol; for serum/plasma sample kits, we provide an optimized serum matrix to mimic the environment of the sample, normalizing assay performance

Reproducible Results You Can Trust

Stability ●●

- All kits are rigorously tested for shipping stability to determine impact on each analyte
- Sample stability is examined at a range of temperatures, and by freeze/thaw

Precision ●●

- Precision of control standards are typically within 15% CV for intra-assay and 20% CV for inter-assay

Accuracy ●

- Spike recovery in sample is performed for each analyte at three concentrations; accuracy is normally within $\pm 30\%$

Dilution Linearity/Parallelism ●

- Diluted sample results must be directly proportional to the concentration of analyte in the sample, typically within $\pm 30\%$

Lot-to-Lot Consistency ●

- A “gold standard” calibrator is defined during assay development and used as a reference for all future production lots

We'll help you choose the right platform and assays to meet your needs today, and we'll partner with you to meet those needs in the future — without compromising quality.

- Screening Assays
- Qualified and Verified Assays

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

MilliporeSigma
400 Summit Drive
Burlington, MA 01803

SigmaAldrich.com



Scan to download
a digital copy.



Keep up to date!

Don't miss out on updates by
signing up for our newsletter:
SigmaAldrich.com/analyteupdate

To place an order or receive technical assistance

Order/Customer Service: **SigmaAldrich.com/order**

Technical Service: **SigmaAldrich.com/techservice**

Safety-related Information: **SigmaAldrich.com/safetycenter**

We have built a unique collection of life science brands with
unrivalled experience in supporting your scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, Millipore, Sigma-Aldrich, Supelco, Milli-Q, SAFC, BioReliance, AAW, Belysa, Conferma, Duolink, FemtoQuest, MAPmate, and MILLIPLEX are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MS_PG7752EN Ver. 9.0
70718
09/2026