

Liquid Solutions MyWay - Europe

Our portfolio of sterile filtered liquids essential for your processes – in stock and ready to ship


Keep your process on the right track with Liquid Solutions MyWay – a portfolio of high-quality sterile filtered liquids in pre-filled single-use bags that are in stock and ready-to-ship. It's that simple: You order, we ship.

The predefined portfolio includes commonly used solutions for your bioprocessing needs. Not only are they in stock, but they have a 12 month shelf life from the date of production. We are adding longer stability studies which will in the future give us a shelf life of 24 months.



Liquid Solutions MyWay: In stock and ready to ship portfolio of high quality sterile filtered liquids

Fast	Convenient	No Compromise on Quality
<ul style="list-style-type: none"> • In stock and ready to ship • Shorter lead times compared to fully custom made solutions 	<ul style="list-style-type: none"> • Order directly from our website • Bags come with standard Colder MPC/MPX connectors 	<ul style="list-style-type: none"> • Ready-to-use bags to mitigate risk and reduce footprint • 12 month shelf life from the date of production (target: 24 months) • EMPROVE® raw materials • Produced in an ISO 9001:2015 certified facility



Robust Packaging Ensures Safe Delivery

Liquid Solutions MyWay single-use bags are ready-to-use and available in sizes from 1 L to 190 L/200 L. All sizes are dual sourced to ensure better availability and supply chain robustness. To safeguard delivery, all shipping containers are ship-tested; hazardous products are packaged in stainless steel while non-hazardous solutions arrive in PET drums.

Exacting Standards Ensure the Highest Quality

All components meet Ph. Eur. requirements, with most of the raw materials also meeting USP. Certificates of analysis (CoA) and safety data sheets (SDS) for the finished goods are available and our manufacturing sites produce in an ISO 9001:2015 certified facility.

Ordering is Fast and Easy

For fast and easy online ordering, visit
SigmaAldrich.com/liquid-solutions-my-way

Liquid Solutions MyWay Product Portfolio

Art. No.	Article	Bag Sizes	Typical Application
55120C	Phosphate Buffered Saline (PBS)	1 L, 20 L, 200 L	Buffer
58300C	Sodium hydroxide solution 0.1 M	1 L, 20 L, 190 L	CIP
58310C	Sodium hydroxide solution 0.5 M	1 L, 20 L, 190 L	CIP
58320C	Sodium hydroxide solution 1 M	1 L, 20 L, 190 L	CIP
58340C	Ethanol solution 20% (v/v)	1 L, 20 L, 200 L	Storage Solution
58345C	Ethanol 20% (v/v), sodium chloride 150 mM	1 L, 20 L, 200 L	Storage Solution

To place an order or receive technical assistance:

Order/Customer Service: SigmaAldrich.com/order

Technical Service: SigmaAldrich.com/techservice

Safety-related Information: SigmaAldrich.com/safetycenter

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

SigmaAldrich.com

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, SAFC, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

