

Calbiochem®

Kit SourceBook™ – 2nd Edition

First choice for Assays & Kits



► *Over 300 Assays & Kits*

► *Over 50 ELISAs*

Categories

- Apoptosis Activity & Detection Kits
- Bioactive Lipid Assays
- Cell Adhesion/Migration Assays
- Cell Cycle–Cell Proliferation Assays
- Cytokine Assays
- Gene & Protein Analysis
- Glycoproteins/Glycobiology
- Growth Factor, Hormone & Receptor Kits
- Heat Shock Protein Assays
- Nitric Oxide/Nitric Oxide Synthase Assays
- Nucleotide Assays
- Oxidative Stress Research Tools
- Phosphatase Assays
- Phospho–Amino Acid Detection Kits
- Protease Assays
- Protein Kinase Assays
- Protein Modification & Labeling Tools
- Protein Quantification Assays
- InhibitorSelect™ 96-Well Protein Kinase Inhibitor Libraries
- WideScreen™ Assays for Luminex™ xMAP® Platform

www.FindYourKit.eu



Cat. No.	Description	No. of Tests	Size	Price £
Apoptosis Activity & Detection Kits				
PF036	Annexin V-Biotin Apoptosis Detection Kit	100	1 ea	229
CBA058	Annexin V-Biotin Apoptosis Detection Kit II	20	1 Kit	107
PF032	Annexin V-FITC Apoptosis Detection Kit	100	1 ea	229
CBA059	Annexin V-FITC Apoptosis Detection Kit II	20	1 Kit	107
CBA060	Annexin V-PE Apoptosis Detection Kit	25	25 Test	107
CBA060	Annexin V-PE Apoptosis Detection Kit	100	100 Test	181
CBA040	Apo-BrdU™ Kit	60	1 Kit	280
CBA041	Apo-Direct™ Kit	50	1 Kit	280
ELISA QIA23	Bcl-2 ELISA Kit	96	1 ea	308
HCS30	BrdU Immunohistochemistry System	50	1 ea	305
218790	Caspase-1 Assay Kit, Colorimetric	100	1 Kit	303
218791	Caspase-1 Assay Kit, Fluorometric	100	1 Kit	303
218792	Caspase-2 Assay Kit, Colorimetric	100	1 Kit	303
218793	Caspase-2 Assay Kit, Fluorometric	96	1 Kit	303
QIA70	Caspase-3 Activity Assay	96	1 Kit	328
HTS02	Caspase-3 Activity Assay for HTS	100	100 Test	263
HTS02	Caspase-3 Activity Assay for HTS	500	500 Test	836
QIA107	Caspase-3 Immunoassay/Activity Kit, Fluorometric	96	1 Kit	314
235419	Caspase-3 Cellular Activity Assay Kit	96	1 Kit	271
235430	Caspase-3 Intracellular Activity Assay Kit I	30	1 Kit	308
235432	Caspase-3 Intracellular Activity Assay Kit II	30	1 Kit	308
QIA91	Caspase-3 Detection Kit (FITC-DEVD-FMK)	100	1 Kit	329
QIA93	Caspase-3 Detection Kit (Red-DEVD-FMK)	100	1 Kit	329
ELISA CBA045	Active Caspase-3 ELISA Kit	96	1 Kit	349
235418	Caspase-3 Inhibitor Screening Assay Kit	96	1 Kit	308
218804	Caspase-5 Assay Kit, Colorimetric	100	1 Kit	303
218802	Caspase-6 Assay Kit, Colorimetric	100	1 Kit	303
QIA108	Caspase-7 Immunoassay/Activity Kit, Fluorometric	96	1 Kit	314
218770	Caspase-8 Assay Kit	96	1 Kit	303
QIA71	Caspase-8 Activity Assay	96	1 Kit	295
QIA113	Caspase-8 Detection Kit (FITC-IETD-FMK)	100	1 Kit	302
QIA114	Caspase-8 Detection Kit (Red-IETD-FMK)	100	1 Kit	302
218824	Caspase-9 Assay Kit, Colorimetric	100	1 Kit	303
QIA72	Caspase-9 Activity Assay	96	1 Kit	328
HTS04	Caspase-9 Activity Assay for HTS	96	100 Test	263
CBA047	Active Caspase-9 Assay Kit	480	1 Kit	275
QIA115	Caspase-9 Detection Kit (FITC-LEHD-FMK)	100	1 Kit	302
QIA116	Caspase-9 Detection Kit (Red-LEHD-FMK)	100	1 Kit	302
218811	Caspase-10 Assay Kit, Fluorometric	96	1 Kit	303
QIA90	Caspase Detection Kit (FITC-VAD-FMK)	100	1 Kit	316
QIA92	Caspase Detection Kit (Red-VAD-FMK)	100	1 Kit	316
ELISA QIA20	Cell Death Detection (Nuclear Matrix Protein) ELISA Kit	96	1 ea	314
HCS17	c-fos IHC System		1 ea	315
ELISA QIA74	Cytochrome c ELISA Kit	96	1 ea	314
QIA87	Cytochrome c Release Apoptosis Assay Kit	100	1 Kit	314
CBA077	InnoCyte™ Cytochrome c Release Kit	50	1 Kit	162
ELISA QIA26	Fas Ligand ELISA Kit	96	1 Kit	316
ELISA QIA26	Fas/APO-1 ELISA Kit	96	1 Kit	238

Calbiochem® – First Choice for Assays & Kits

Cat. No.	Description	No. of Tests	Size	Price £
QIA33	TdT FragEL™ DNA Fragmentation Detection Kit	50	1 EA	302
QIA39	Fluorescein FragEL™ DNA Fragmentation Detection Kit	50	1 EA	261
QIA21	Klenow FragEL™ DNA Fragmentation Detection Kit	50	1 EA	302
QIA89	Glutathione Apoptosis Detection Kit, Fluorometric	100	1 Kit	274
CBA014	InnoZyme™ Granzyme B Immunocapture Activity Assay Kit	96	1 Kit	313
475866	MitoCapture™ Apoptosis Detection Kit	100	1 Kit	254
 QIA25	Nucleosome ELISA Kit	96	1 EA	210
 QIA03	Mutant p53 ELISA Kit	36 or 42	1 EA	306
 QIA07	p53 ELISA Kit (Pantropic)	96	1 EA	306
 QIA53	p53 ELISA ^{PLUS} (Autoantibody) Kit	96	1 EA	350
 CBA009	PhosphoDetect™ p53 (pSer ¹⁵) ELISA Kit	96	1 Kit	369
512729	PARP Cleavage Detection Kit		1 Kit	171
AM41	Suicide Track™ DNA Ladder Isolation Kit	25	1 EA	145
 CBA048	Survivin ELISA Kit	182	1 Kit	263
 CBA049	XIAP ELISA Kit	182	1 Kit	233
Bioactive Lipid Assays				
428901	Cholesterol/Cholesteryl Ester Quantitation Kit	100	1 Kit	210
428900	LCAT Activity Assay Kit, Fluorometric	100	1 Kit	308
Cell Adhesion/Migration Assays				
QIA129	Innocyte™ Cell Invasion Assay Kit (24-well)	12	1 Kit	281
 QIA112	Galectin-3 ELISA Kit	96	1 Kit	348
CBA037	Innocyte™ α2, α3, α5 Integrin Detection Kit	96	1 Kit	204
CBA017	Innocyte™ Cell Migration Assay, 24-well plate	12	1 Kit	250
CBA010	InnoCyte™ Cell Migration Assay, 96-well	96	1 Kit	255
CBA011	InnoCyte™ ECM Cell Adhesion Assay, Fibronectin	72	1 Kit	126
CBA012	InnoCyte™ ECM Cell Adhesion Assay, Vitronectin	72	1 Kit	184
CBA013	InnoCyte™ ECM Cell Adhesion Assay, Collagen Type IV	72	1 Kit	151
CBA025	InnoCyte™ ECM Cell Adhesion Assay, Laminin/BMC	72	1 Kit	193
CBA026	InnoCyte™ ECM Cell Adhesion Assay, Galectin-1/-3	72	1 Kit	204
CBA043	InnoCyte™ Laminin-based 96-Well Cell Invasion Assay	96	1 Kit	247
CBA083	InnoCyte™ ECM Cell Adhesion Assay	96	1 Kit	273
Cell Cycle-Cell Proliferation Assays				
QIA58	BrdU Cell Proliferation Assay	200	200 Test	133
QIA58	BrdU Cell Proliferation Assay	1000	1000 Test	312
HTS01	BrdU Cell Proliferation Assay, HTS	200	200 Test	142
HTS01	BrdU Cell Proliferation Assay, HTS	1000	1000 Test	321
QIA76	Live/Dead Double Staining Kit	100	100 Test	186
 QIA59	PCNA ELISA Kit	96	1 EA	306
QIA127	Rapid Cell Proliferation Kit	500	1 Kit	129
QIA117	Senescence Detection Kit	250	1 Kit	183
QIA128	Ultrasensitive Cell Proliferation Assay Kit, Fluorogenic	500	1 Kit	143
 QIA18	p21 ^{WAF1} ELISA Kit	96	1 EA	411
Cytokine Assays				
QIA68	IL-8 Accucyte® EIA Kit, Human	96	1 EA	244
Gene & Protein Analysis				
71377-3	NoShift™ Transcription Factor Assay Kit	96	1 Kit	290
71378-3	NoShift™ Transcription Factor Assay Kit Plus NucBuster™	96	1 Kit	377
71518-3	NoShift™ NF-κB (P65) Reagents	100	1 Kit	165

NEW

Calbiochem® - First Choice for Assays & Kits

Cat. No.	Description	No. of Tests	Size	Price £
71678-3	NoShift™ II EGR Transcription Factor Assay	96	1 Kit	350
71679-3	NoShift™ II NF-1 Transcription Factor Assay	96	1 Kit	350
71680-3	NoShift™ II NF-κB Transcription Factor Assay	96	1 Kit	350
71681-3	NoShift™ II PPAR Transcription Factor Assay	96	1 Kit	350
Glycoproteins/Glycobiology				
362280	Glycoprotein Deglycosylation Kit	10	1 Kit	334
Growth Factor, Hormone & Receptor Kits				
 CBA018	Epidermal Growth Factor Receptor (EGF-R) ELISA Kit	96	1 Kit	369
 QIA64	Basic FGF ELISA Kit	96	1 ea	416
 QIA61	TGF-α ELISA Kit	96	1 ea	306
 QIA51	VEGF ELISA Kit, Human	100	1 ea	367
 QIA52	VEGF ELISA Kit, Mouse	192	1 ea	367
Heat Shock Protein Assays				
 QIA119	Hsp27 ELISA Kit	96	1 Kit	352
 CBA078	Phosphodetect™ HSP27 (pSer ^{78/82}) ELISA Kit	96	1 Kit	252
 CBA057	Hsp90α ELISA Kit	96	1 Kit	320
Nitric Oxide/Nitric Oxide Synthase Assays				
482650	Nitric Oxide Assay Kit, Colorimetric	80	1 Kit	228
482655	Nitric Oxide Assay Kit, Fluorometric	192-288	1 Kit	263
482700	Nitric Oxide Synthase Assay Kit	50	1 Kit	310
482702	Nitric Oxide Synthase Assay Kit, Colorimetric	96	1 Kit	263
Nucleotide Assays				
119107	ATP Assay Kit	200	1 Kit	225
116811	cAMP Direct Immunoassay Kit, Colorimetric	100	1 Kit	308
Oxidative Stress Research Tools				
219265	Catalase Assay Kit	96	1 Kit	259
 CBA053	COX-2 ELISA Kit	96	1 Kit	340
 QIA97	Cu/Zn Superoxide Dismutase ELISA	96	1 Kit	314
353919	Glutathione Peroxidase Assay Kit	96	1 Kit	259
354102	Glutathione Assay Kit	100	100 Test	313
354103	Glutathione Assay Kit II	96	1 Kit	241
359963	Glutathione Reductase Assay Kit II	96	1 Kit	241
371757	GSH/GSSG Ratio Assay Kit	200	1 Kit	472
 427275	Lactoferrin ELISA Kit	96	1 Kit	366
437639	Lipid Hydroperoxide (LPO) Assay Kit	100	1 Kit	243
CBA024	InnoZyme™ Myeloperoxidase Activity Kit	96	1 Kit	335
 475919	Myeloperoxidase ELISA Kit	96	1 Kit	384
500095	OxyDNA Assay Kit, Fluorometric	50	1 Kit	451
574590	Superoxide Anion Detection Kit, Chemiluminescent	100	1 Kit	306
574601	Superoxide Dismutase Assay Kit II	96	1 Kit	241
615700	Total Antioxidant Status Assay Kit	50	1 Kit	592
Phosphatase Assays				
207005	Calcineurin Assay Kit, Colorimetric	100	1 Kit	314
207007	Calcineurin Cellular Activity Assay Kit, Colorimetric	192	1 Kit	272
539736	Protein Tyrosine Phosphatase 1B Assay Kit, Colorimetric	96	1 Kit	276
539742	PTP Assay Kit, Fluorometric	96	1 Kit	204
Phospho-Amino Acid Detection Kits				
525282	Phosphoserine Detection Kit		1 Kit	310

NEW

Calbiochem® – First Choice for Assays & Kits

Cat. No.	Description	No. of Tests	Size	Price £
525288	Phosphothreonine Detection Kit		1 Kit	263
525322	Phosphotyrosine Detection Kit		1 Kit	241
Protease Assays				
565785	β-Secretase Activity Assay Kit, Fluorogenic	100	1 Kit	306
CBA054	InnoZyme™ Calpain1/2 Activity Assay Kit	96	1 Kit	227
QIA120	Calpain Activity Assay Kit, Fluorogenic	96	1 Kit	299
CBA003	InnoZyme™ Gelatinase Activity Assay Kit, Fluorogenic	96	1 Kit	324
CBA001	InnoZyme™ Cathepsin B Activity Assay Kit, Fluorogenic	96	1 Kit	229
219346	Cathepsin B Detection Kit		1 Kit	228
CBA002	InnoZyme™ Cathepsin D Immunocapture Activity Assay Kit	96	1 Kit	283
ELISA QIA29	Cathepsin D, Rapid Format ELISA Kit	96	1 EA	306
219347	Cathepsin K Detection Kit		1 Kit	228
CBA023	InnoZyme™ Cathepsin L Activity Kit, Fluorogenic	96	1 Kit	310
ELISA QIA94	Cathepsin L ELISA Kit	96	1 Kit	316
CBA050	InnoZyme™ Cathepsin S Activity Assay Kit	96	1 Kit	275
ELISA QIA55	MMP-1 ELISA	96	1 ea	297
ELISA QIA63	MMP-2 ELISA Kit	34 or 41	1 ea	416
ELISA QIA73	MMP-3 ELISA Kit	34 or 41	1 ea	362
ELISA QIA56	MMP-9 ELISA Kit	34 or 41	1 ea	411
ELISA QIA130	MMP-13 ELISA Kit	96	1 Kit	338
539125	Protease Assay Kit	200	1 Kit	184
539158	20S Proteasome Assay Kit, SDS-Activated	30 or 150	1 Kit	271
539159	26S Proteasome Degradation Kit		1 Kit	262
539176	Proteasome Isolation Kit, Human		1 Kit	262
ELISA QIA54	TIMP-1 ELISA Kit	96	1 ea	297
ELISA QIA40	TIMP-2 ELISA Kit	96	1 ea	389
ELISA QIA96	PMN Elastase/α1-PI Complex ELISA Kit	96	1 Kit	306
CBA016	InnoZyme™ Human Neutrophil Elastase Immunocapture Activity Kit	96	1 Kit	313
ELISA QIA99	Soluble t-PA ELISA Kit	96	1 Kit	308
CBA042	InnoZyme™ TACE Activity Kit	96	1 Kit	274
662096	Ubiquitin Protein Conjugation Kit		1 Kit	262
662200	Ubiquitinated Protein Enrichment Kit		1 Kit	262
QIA125	Urokinase Activity Assay Kit, Fluorogenic	100	1 Kit	226
Protein Kinase Assays				
ELISA CBA019	K-LISA™ Akt Activity Kit	96	1 Kit	263
124007	Akt Activity Immunoassay Kit	40	1 Kit	286
ELISA CBA004	PhosphoDetect™ Akt (pThr ³⁰⁸) ELISA Kit	96	1 Kit	369
ELISA CBA005	PhosphoDetect™ Akt (pSer ⁴⁷³) ELISA Kit	96	1 Kit	369
ELISA CBA036	PhosphoDetect™ ATF2 (pThr ^{69/71}) ELISA Kit	96	1 Kit	369
ELISA CBA051	K-LISA™ Aurora Kinase Activity Kit	96	1 Kit	237
CBA052	IP-Activity™ Bcr-Abl/c-Abl Kit		1 Kit	98
ELISA CBA020	K-LISA™ Checkpoint Activity Kit	96	1 Kit	263
ELISA QIA10	c-erbB2/c-neu Rapid Format ELISA Kit	96	1 ea	416
ELISA CBA074	c-Kit ELISA Kit	96	1 Kit	373
ELISA CBA075	PhosphoDetect™ c-Kit (pTyr ⁸²³) ELISA Kit	96	1 Kit	373
ELISA CBA073	PhosphoDetect™ c-Met (pTyr ^{1230/1234/1235}) ELISA Kit	96	1 Kit	373



Calbiochem® – First Choice for Assays & Kits

	Cat. No.	Description	No. of Tests	Size	Price £
ELISA	CBA071	CREB ELISA Kit	96	1 Kit	373
ELISA	CBA072	PhosphoDetect™ CREB (pSer ¹³³) ELISA Kit	96	1 Kit	373
ELISA	CBA032	ERK 1/2 ELISA Kit	96	1 Kit	369
ELISA	CBA006	PhosphoDetect™ ERK1/2 (pThr ¹⁸⁵ /pTyr ¹⁸⁷) ELISA Kit	96	1 Kit	369
ELISA	CBA063	FAK ELISA Kit	96	1 Kit	373
ELISA	CBA062	PhosphoDetect™ FAK (pTyr ³⁹⁷) ELISA Kit	96	1 Kit	373
ELISA	CBA068	GSK-3β ELISA Kit	96	1 Kit	373
ELISA	CBA069	PhosphoDetect™ GSK-3β (pSer ⁹) ELISA Kit	96	1 Kit	373
ELISA	CBA064	IκBα ELISA Kit	96	1 Kit	373
ELISA	CBA065	PhosphoDetect™ IκBα (pSer ³²) ELISA Kit	96	1 Kit	373
	CBA044	K-LISA™ IKKβ Inhibitor Screening Kit	96	1 Kit	237
ELISA	CBA039	Insulin Receptor (β-subunit) ELISA Kit	96	1 Kit	369
ELISA	CBA038	PhosphoDetect™ Insulin Receptor (pTyr ^{1162/1163}) ELISA Kit	96	1 Kit	369
	CBA079	InnoZyme™ Insulysin Immunocapture Activity Kit	96	1 Kit	275
ELISA	CBA033	JNK 1/2 ELISA Kit	96	1 Kit	369
ELISA	CBA007	PhosphoDetect™ JNK1/2 (pThr ¹⁸³ /pTyr ¹⁸⁵) ELISA Kit	96	1 Kit	369
	420115	JNK Activity Immunoassay Kit	40	1 Kit	313
ELISA	CBA070	MEK1 ELISA Kit	96	1 Kit	373
ELISA	CBA030	PhosphoDetect™ MEK1 (pSer ^{218/222}) ELISA Kit	96	1 Kit	369
ELISA	CBA055	K-LISA™ mTOR Activity Kit	96	1 Kit	255
ELISA	CBA029	p38 MAP Kinase ELISA Kit	96	1 Kit	369
ELISA	CBA008	PhosphoDetect™ p38 MAP Kinase (pThr ¹⁸⁰ /pTyr ¹⁸²) ELISA Kit	96	1 Kit	369
	539484	Protein Kinase Assay Kit, Non-Radioactive	96	1 Kit	369
ELISA	CBA066	PRAS40 ELISA Kit	96	1 Kit	373
ELISA	CBA067	PhosphoDetect™ PRAS40 (pThr ²⁴⁶) ELISA Kit	96	1 Kit	373
ELISA	CBA027	Src ELISA Kit	96	1 Kit	369
	CBA056	IP-Activity™ Src Kit		1 Kit	98
ELISA	CBA028	PhosphoDetect™ Src (pTyr ⁴¹⁸) ELISA Kit	96	1 Kit	369
ELISA	CBA034	STAT1 ELISA Kit	96	1 Kit	369
ELISA	CBA035	PhosphoDetect™ STAT1 (pTyr ⁷⁰¹) ELISA Kit	96	1 Kit	369
Protein Modification & Labeling Tools					
	343210	FITC Labeling Kit	5	1 Kit	173
	XHC01	Mouse UniTect™ ABC Kit	800	1 Kit	119
	XHC02	Rabbit UniTect™ ABC Kit	800	1 Kit	119
Protein Quantification Assays					
	219468	CB-Protein Assay™ Kit	500 or 5000	1 Kit	111
	488250	Non-Interfering Protein Assay™ Kit	500	1 Kit	190
Protein Sample Preparation Assays					
	122640	ProteoExtract® Albumin Removal Kit	12	1 Kit	142
	122641	ProteoExtract® Albumin Removal Kit, Maxi	20	1 Kit	148
	122642	ProteoExtract® Albumin/IgG Removal Kit	12	1 Kit	152
	122643	ProteoExtract® Albumin/IgG Removal Kit, Maxi	20	1 Kit	169
	444810	ProteoExtract® Native Membrane Protein Extraction Kit	20	1 Kit	238
	539180	ProteoExtract® Protein Precipitation Kit	200	1 Kit	161
	539182	ProteoExtract® Tryptic Cleavage Modification Kit	100	1 Kit	83
	539779	ProteoExtract® Complete Mammalian Proteome Extraction Kit	20	1 Kit	208

Calbiochem® – First Choice for Assays & Kits

Cat. No.	Description	No. of Tests	Size	Price £	
539779	ProteoExtract® Complete Mammalian Proteome Extraction Kit	20	1 Kit	208	
539789	ProteoExtract® Partial Mammalian Proteome Extraction Kit	20	1 Kit	242	
539790	ProteoExtract® Subcellular Proteome Extraction Kit	20	1 Kit	266	
539791	ProteoExtract® Subcellular Proteome Extraction Kit, Mini	10	1 Kit	141	
650212	ProteoExtract® All-in-One Trypsin Digestion Kit	100	1 Kit	186	
QIA88	ProteoExtract® Cytosol/Mitochondria Fractionation Kit	100	1 Kit	227	
539720	ProteoExtract® Tissue Dissociation Buffer Kit	10	1 Kit	50	
539723	ProteoExtract® Phosphopeptide Enrichment SCIMAC Kit	100	1 Kit	144	
539722	ProteoExtract® Phosphopeptide Enrichment TiO ₂ Kit	100	1 Kit	149	
71772-3	ProteoExtract® Transmembrane Protein Extraction Kit	20	1 Kit	173	
InhibitorSelect™ Protein Kinase Inhibitor Libraries					
539744	InhibitorSelect™ 96-Well Protein Kinase Inhibitor Library I		1 ea	2195	
539745	InhibitorSelect™ 96-Well Protein Kinase Inhibitor Library II		1 ea	2195	
WideScreen™ Assays for Luminex® xMAP® Platform					
Widescreen™ RTK Assays					
71942-3	WideScreen™ RTK Total Assay Complete Kit	96	1 Kit	1295	
71943-3	WideScreen™ RTK pTyr Assay Complete Kit	96	1 Kit	1250	
71783-3	WideScreen™ Reagent Kit		1 Kit	120	
71926-3	Cell Extraction Kit		1 Kit	48	
71924-3	RTK Total 7-Plex	100	1 Kit	1170	
71925-3	RTK pTyr 7-Plex	100	1 Kit	1105	
71728-3	IGF-1R Total Bead Kit	100	1 Kit	180	
71929-3	EGFR Total Bead Kit	100	1 Kit	180	
71930-3	HGFR Total Bead Kit	100	1 Kit	180	
71931-3	PDGFRβ Total Bead Kit	100	1 Kit	180	
71932-3	HER-2 Total Bead Kit	100	1 Kit	180	
71933-3	VEGFR2 Total Bead Kit	100	1 Kit	180	
71934-3	Tie-2 Total Bead Kit	100	1 Kit	180	
71935-3	EGFR pTyr Bead Kit	100	1 Kit	170	
71936-3	IGF-1R pTyr Bead Kit	100	1 Kit	170	
71937-3	HGFR pTyr Bead Kit	100	1 Kit	170	
71938-3	PDGFRβ pTyr Bead Kit	100	1 Kit	170	
71939-3	HER-2 pTyr Bead Kit	100	1 Kit	170	
71940-3	VEGFR2 pTyr Bead Kit	100	1 Kit	170	
71941-3	Tie-2 pTyr Bead Kit	100	1 Kit	170	
WideScreen™ EpiTag™ Assays					
71782-3	WideScreen™ EpiTag™ ERK Pathway Panel I Complete Kit	96	1 Kit	1195	
71891-3	WideScreen™ EpiTag™ ERK Pathway Panel II Complete Kit	96	1 Kit	1195	
71784-3	EpiTag™ Cell Lysis Kit		1 Kit	38	
71785-3	EpiTag™ Digest Kit		1 Kit	38	
71791-3	EpiTag™ ERK Pathway Standards Mix		1 Kit	30	
71794-3	EpiTag™ Combo Reagent Kit		1 Kit	180	
71726-3	EpiTag™ ERK Panel I	100	1 Kit	995	
71890-3	EpiTag™ ERK Panel II	100	1 Kit	995	
71829-3	EpiTag™ ERK1 Total Bead Kit	100	1 Kit	185	
71830-3	EpiTag™ ERK2 Total Bead Kit	100	1 Kit	185	
71831-3	EpiTag™ B-Raf (pS446) Bead Kit	100	1 Kit	185	
71832-3	EpiTag™ MEK1/2 (pS217, pS221/ pS222, pS226) Bead Kit	100	1 Kit	185	

Cat. No.	Description	No. of Tests	Size	Price £	
71834-3	EpiTag™ ERK2 (pT185, pY187) Bead Kit	100	1 Kit	185	
71835-3	EpiTag™ MEK1 Total Bead Kit	100	1 Kit	185	
71836-3	EpiTag™ MEK2 Total Bead Kit	100	1 Kit	185	
71837-3	EpiTag™ STAT1 Total Bead Kit	100	1 Kit	185	
71838-3	EpiTag™ Raf1 Total Bead Kit	100	1 Kit	185	
71839-3	EpiTag™ Raf1 (pS338) Bead Kit	100	1 Kit	185	

Online Resources

www.FindYourKit.eu

- Browse kits by category
- Browse kits by type
- Kit comparison tables
- Technical information

www.FindYourInhibitor.eu

- Browse inhibitors by category
- Download substructure searchable database
- Technical tips

www.FindYourAntibody.eu

- Search for antibodies by target antigen
- Browse antibodies by name
- Browse antibodies by category

www.FindYourPathway.eu

- More than seventeen pathways available with new pathways continually added
- Detailed and up-to-date schematics

For more information or to order
Calbiochem products, contact

Merck Chemicals Ltd.

Boulevard Industrial Park
Padge Road, Beeston
Nottingham NG9 2JR

Visit our website

www.merckbio.eu

Freephone

Fax

E-mail

Orders / Enquiries

0800 622 935

0115 943 0951

customer.service@merckbio.eu

Technical Support

Freephone

Fax

E-mail

0800 622 935

0115 943 0951

techservice@merckbio.eu